

## PRODUCT FEATURES

### 1. MORE CONVENIENT!

- Completely portable, easy and convenient to use
- You don't have to hold the device (unlike a laser hair comb) and completely fits in almost any hat
- Can be used anywhere and while engaging in other daily activities, such as at home, walking, driving...anywhere

### 2. HIGHLY EFFECTIVE!

- Can be used in combination with other therapies, such as Minoxidil, and hair transplant
- More laser diodes and higher power output than other portable home use product
- 272 individual 650nm low level laser diodes, no LED light used

### 3. CHEAPER!

- More cost-effective than other hair growth devices applying LLLT
- You don't need to spend much visit a salon or physician's office multiple times per week to receive treatment.
- Affordable

## COMPLETE STANDARD PACKAGE

### Laser Therapy Cap includes

- 1 laser cap
- 1 battery pack
- 1 wall charger
- 1 baseball hat
- 1 carrying case
- 1 user instruction



Introducing more and better medical cosmetic and dermatological treatment equipments is our pursuit and ideal of a lifetime  
We will accept the challenge bravely!



## HairPro

LOW LEVEL LASER THERAPY  
FOR THINNING HAIR



## TREATMENT PRINCIPLE

1.The device utilizes low level laser therapy which can enhance the activity of Adenosine Triphosphate (ATP) in subcutaneous tissue and accelerate blood circulation of the scalp, and then the oxygen and nutrient can be sent to hair follicles promptly, as a result, the efficacy of halting the progression hair loss can be achieved.

2.Low Level laser can dredge sebaceous gland ducts to decrease and control oil secretion so as to improve hair condition.

3.Low level laser therapy can effectively stimulate most hair follicles which are in the state of hibernation to active status and improve hair loss in telogen and too short hair growth period, to achieve the effect of reducing hair loss and hair re-growth.

## INDICATIONS

Hair growth  
Hair restoration  
Anti-hair loss



## LASER THERAPY CAP

It is a portable dome-shaped device that consists of individual 650nm laser diodes, each with 5mW of power. Automatic session programmed for 30-minute sessions, with commencement and end indicated by distinct audible sounds, the cap can be used virtually anywhere—at home, while travelling, driving, watching TV, or even on a jog. It is worn under any hat of your preference for half an hour every other day or as otherwise prescribed by your physician.

### LASER THERAPY RESULTS



Powered by a rechargeable battery pack, the HairPro battery pack can be slipped into a pocket and offers patients the low level laser power at home.

NOTICE: The HairPro is a cosmetic device for personal home use.



## MOBILE LASER THERAPY CAP MODEL

### 272 Laser Hair Cap

Laser Wavelength: 650nm ± 5nm  
Number of lasers: 272  
Diode output: 5mW  
Total output: 1360mW



### 81 Laser Hair Cap

Laser Wavelength: 650nm ± 5nm  
Number of lasers: 81  
Diode output: 5mW  
Total output: 405mW

