




MICRO ENDOSCOPES

MEDICAL GRADE MICRO ENDOSCOPE



- 0.8 mm DIAMETER, 75 cm LENGTH FLEXIBLE SCOPE
- IMAGE TRANSMISSION: 6,000 PIXELS
- DIRECTION OF VIEW: 0°
- FIELD OF VIEW : 70°
- INTEGRATED FIBER OPTIC ILLUMINATION
- STANDARD 32 MM DIA ENDOSCOPE EYEPIECE

VIDEO MICRO ENDOSCOPE




- 0.6 mm DIAMETER, 75 cm LENGTH FLEXIBLE SCOPE
- IMAGE TRANSMISSION: 6,000 PIXELS
- DIRECTION OF VIEW: 0°
- FIELD OF VIEW : 70°
- INTEGRATED FIBER OPTIC ILLUMINATION
- QUICK SICONNECT VIDEO CAMERA MOUNT

VIDEO CAMERA

- 1/4" COLOR CCD CAMERA
- BUILT-IN WHITE ADJUSTABLE LED ILLUMINATION
- WHITE BALANCE BUTTON
- BUILT-IN VIDEO COUPLER WITH FOCUSING
- COMPACT ERGONOMIC DESIGN
- NTSC OR USB OOUTPUT

VIDEO SEMI-RIGID MICRO ENDOSCOPE



- 0.9 mm DIAMETER, 150 mm LONG SEMI-RIGID SCOPE
- BUILT-IN FIBER OPTIC ILLUMINATION
- ONE WORKING CHANNEL

VIDEO CAMERA OPTIONS

- **1/4" COLOR CCD REMOTE HEAD CAMERA**
- CAMERA CONTROL UNIT WITH WHITE BALANCE CONTROL
- SEALED CAMERA HOUSING
- VIDEO COUPLER WITH FOCUSING
- NTSC OUPUT



- **1/6" COLOR CCD REMOTE HEAD CAMERA**
- CAMERA CONTROL UNIT WITH WHITE BALANCE CONTROL
- BUILT-IN FIBER OPTIC ILLUMINATION
- REMOVABLE CAMERA HEAD WITH OPTICS
- NTSC OUTPUT