



**AGS-Medical**



## Breathing Circuit For Anaesthesia

### Coaxial Circuit/Bain's Circuit

- ▶ Low flow anaesthesia in mind
- ▶ Offering warm
- ▶ More humid inspiratory gases




SPECIFICATION	MODEL	LENGTH	DESCRIPTION	Packing (carton)
Disposable Coaxial Circuits	Standard	1.5m/1,8m	Bain's circuit, w/ wo breathing bag	30pcs
	With limb	1.5m/1,8m	Bain's circuit with 0.6m limb, w/ wo breathing bag	30pcs

Contact us for more specifications

### Corrugated Circuit/Standard Circuit

- ▶ Good flexibility, bending resistance
- ▶ Consist of corrugated hoses
- ▶ Swivel elbow with the luer port and cap



**Rotatable**

SPECIFICATION	MODEL	LENGTH	DESCRIPTION	Packing (carton)
Disposable Corrugated Anesthesia Circuits	Adult	1.0m	Φ22mm, w/ wo breathing bag	25pcs
		1.5m	Φ22mm, w/ wo breathing bag	20pcs
		1.8m	Φ22mm, w/ wo breathing bag	20pcs
	Paediatric	1.0m	Φ19mm, w/ wo breathing bag	30pcs
		1.5m	Φ19mm, w/ wo breathing bag	30pcs

Contact us for more specifications

**Ph: 505-550-6501 & 505-565-5102, Fx: 505-814-5778**

**Em: sales@agsmedical.com, Web: www.agsmedical.com**



AGS-Medical




## Collapsible Circuit/Telescopic Circuit

- ▶ Adjustable length
- ▶ High transparency and odor free
- ▶ Coaxial circuit




### Collapsible circuit



SPECIFICATION	MODEL	LENGTH	DESCRIPTION	Packing (carton)
Disposable Collapsible Anesthesia Circuits	Adult	1.8m/2.4m	F22mm ,w/ wo breathing bag	50pcs
	Paediatric	1.5m/1.8m	F19mm ,w/ wo breathing bag	50pcs

## Limb Circuit/CPAP



SPECIFICATION	MODEL	DESCRIPTION
Limb circuit	Corrugated circuit	Corrugated extra limb
	Collapsible circuit	Expandable extra limb
	Smooth Bore circuit	Smooth Bore extra limb

Contact us for more specifications

Ph: 505-550-6501 & 505-565-5102, Fx: 505-814-5778

Em: sales@agsmedical.com, Web: www.agsmedical.com





AGS-Medical



## Breathing Circuit For Ventilator Standard Circuit

- ▶ High transparency and odor free
- ▶ Good flexibility, bending resistance
- ▶ Standard connector




SPECIFICATION	MODEL	LENGTH	DESCRIPTION	Packing (carton)
Disposable Corrugated Breathing Circuits	One water trap	1,6m	With one water trap breathing circuit and 0.6m limb	20pcs
		1,8m	With one water trap breathing circuit and 0.6m limb	20pcs
	Two water traps	1,6m	With two water traps breathing circuit and 0.6m limb	20pcs
		1,8m	With two water traps breathing circuit and 0.6m limb	20pcs

Contact us for more specifications

## Smooth Circuit

- ▶ Improve flow characteristics
- ▶ Low tubing compliance



SPECIFICATION	MODEL	LENGTH	DESCRIPTION	Packing (carton)
Disposable Smooth Bore Breathing Circuits	One water trap	1,6m	With one water trap breathing circuit and 0.6m limb	20pcs
		1,8m	With one water trap breathing circuit and 0.6m limb	20pcs
	Two water traps	1,6m	With two water traps breathing circuit and 0.6m limb	20pcs
		1,8m	With two water traps breathing circuit and 0.6m limb	20pcs

Contact us for more specifications

Ph: 505-550-6501 & 505-565-5102, Fx: 505-814-5778


Em: sales@agsmedical.com, Web: www.agsmedical.com



AGS-Medical



## Accessories



Ph: 505-550-6501 & 505-565-5102, Fx: 505-814-5778

Em: [sales@agsmedical.com](mailto:sales@agsmedical.com), Web: [www.agsmedical.com](http://www.agsmedical.com)