




Web: • <http://www.agsmedical.com> • Email: sales@agsmedical.com • Skype: agstech1
 Tel: (505) 550-6501 • Fax: (505) 814-5778


EVERVU VIDEO CAMERA SYSTEM FOR ENDOSCOPIC APPLICATIONS

The system represents the new generation of endoscopic video cameras combining the high efficiency LED illumination with full digital signal processing CCD camera.




Light Source

LED Output	80 W
Light Color Temperature	5600°K
Power Supply Voltage	90-264 V AC, Single Phase
Frequency	50*60 Hz
Power Consumption	175 W
Class	1 Type B
LED Service Life	50,000 hrs (approx.)
Cooling	Forced air
Overheat Protection	Automatic LED Power-OFF
Dimmer	Electronic
Fiber Optic Connector	Built-in Revolving Turret with Multiple Industry Standard Adapters
Dimensions	202(W) x 290(D) x 95(H) mm
Weight	3.0 kg.



Video Camera

Sensor	1/3" Color CCD
Resolution	525 TV lines
Shutter Speed	Up to 1/100,000 sec.
Min. Illumination	2.0 Lux F1.4
Video Output	Composite Video (NTSC), 1Vpp/75 ohm, BNC S-Video (Y-C), Mini-DIN Jack
Handpiece Controls	White Balance Anti-Moire (for flexible scopes) Window Size Photo Capture



Video Coupler

Focusing	Manual
Focal Length	Fixed : 25, 28, 30, 35, 40, 45 mm
Camera Mount	C-Mount
Sealing	Sealed, Soakable
Clamp	Spring Loaded, Fits 32 mm Eyepiece

SEE US ON THE WEB AT: www.agsmedical.com