

## PRODUCT OVERVIEW

### Multi-DOA test profile

PRODUCT	SPECIMEN	TYPE	SENSITIVITY
Morphine	urine	strip/cassette	300ng/ml
Methyl amphetamine	urine	strip/cassette	1000ng/ml
Amphetamine	urine	strip/cassette	1000ng/ml
Cocaine	urine	strip/cassette	300ng/ml
Ketamine	urine	strip/cassette	1000ng/ml
Methodone	urine	strip/cassette	300ng/ml
THC	urine	strip/cassette	50ng/ml
Barbital	urine	strip/cassette	300ng/ml
MDMA	urine	strip/cassette	500ng/ml
Benzodiazepams	urine	strip/cassette	300ng/ml
Bupivacaine	urine	strip/cassette	10ng/ml

### Infectious Disease test profile

PRODUCT	SPECIMEN	TYPE	SENSITIVITY
HIV 1/2	serum/plasma/whole blood	strip/cassette	4 NCU/ml
Syphilis-TP	serum/plasma/whole blood	strip/cassette	2 NCU/ml
Malaria PF	whole blood	strip/cassette	50 parasites/ul
Malaria PV/PF	whole blood	strip/cassette	50 parasites/ul
Malaria Pan/Pf	whole blood	strip/cassette	50 parasites/ul
Tuberculosis	serum/plasma/whole blood	strip/cassette	N/A
Dengue IgG/IgM	serum/plasma/whole blood	strip/cassette	N/A
Dengue NS1	serum/plasma/whole blood	strip/cassette	N/A
Herpes simplex	serum/plasma	strip/cassette	N/A
Mycoplasma	serum/plasma	strip/cassette	N/A
Chlamydia	Cervical Mucus	strip/cassette	N/A
Gonorrhea	Cervical Mucus	strip/cassette	N/A
H.pylori antibody	serum/plasma/whole blood	strip/cassette	N/A
Typhoid	serum/plasma/	cassette	

### Tumor Marker test profile

PRODUCT	SPECIMEN	TYPE	SENSITIVITY
AFP(Alpha-fetoprotein)	serum/plasma/whole blood	strip/cassette	20ng/ml
CEA(Carcinoembryonic antigen)	serum/plasma/whole blood	strip/cassette	5ng/ml
FOB(Fecal occult blood)	Feces	strip/cassette	40ng/ml
PSA(Prostate specific antigen)	serum/plasma/whole blood	strip/cassette	4ng/ml

### Fertility&Eugenic Testing Product Profile

PRODUCT	SPECIMEN	TYPE	SENSITIVITY
Pregnancy(HCG)	Urine/serum	strip/cassette/midstream	25mIU/ml
Ovulation(LH)	Urine	strip/cassette/midstream	20mIU/ml
FSH(Follicle-stimulating Hormone)	urine	strip/cassette	25mIU/ml
fFN(Fetal Fibronectin)	Cervical Mucus	strip/cassette	50ng/ml
HSV 1/HSV 2 antibody	serum/plasma	strip/cassette	N/A
TOXO antibody	serum/plasma	strip/cassette	N/A
Rubella Virus antibody	serum/plasma	strip/cassette	N/A
Cytomegalovirus antibody	serum/plasma	strip/cassette	N/A
TORCH-IgM(COMB-5)	serum/plasma	cassette	N/A

### Liver Disease test profile

PRODUCT	SPECIMEN	TYPE	SENSITIVITY
Hepatitis A IgM	serum/plasma/whole blood	strip/cassette	4 NCU/ml
HBsAg	serum/plasma/whole blood	strip/cassette	0.5ng/ml
HBsAb	serum/plasma/whole blood	strip/cassette	10mIU/ml
HbeAg	serum/plasma/whole blood	strip/cassette	2NCU/ml
HBeAb	serum/plasma/whole blood	strip/cassette	N/A
HBcAb	serum/plasma/whole blood	strip/cassette	2 NCU/ml
Multi-HBV	serum/plasma/whole blood	5 combo	N/A
HCV IgM	serum/plasma/whole blood	strip/cassette	2NCU/ml
HEV IgM	serum/plasma/whole blood	strip/cassette	2NCU/ml
HAV/HEV IgM	serum/plasma/whole blood	2 combo	4 NCU/ml

AGS-Medical, Ph: 505-550-6501 & 505-565-5102, Fx: 505-814-5778, Em: [sales@agsmedical.com](mailto:sales@agsmedical.com), Web: <http://www.agsmedical.com>

Reference Code: OICASSHANDONGHARMOWELLOICAS