



AGS-Medical



## Anesthesia Mask

### Type C

- ▶ Skid resistance, easy handling
- ▶ Soft and flexible–transparent air cushion
- ▶ Standard connector
- ▶ Removable hook ring



| SPECIFICATION            | MODEL            | CODE NUMBER | DESCRIPTION                         | Packing (carton) |
|--------------------------|------------------|-------------|-------------------------------------|------------------|
| Anesthesia Mask (Type C) | Neonate (1#)     | 5401310     | 15mm OD,skid resistance, inflatable | 30pcs            |
|                          | Infant (2#)      | 5401320     | 15mm OD,skid resistance,inflatable  | 30pcs            |
|                          | Child (3#)       | 5401330     | 22mm OD,skid resistance, inflatable | 70pcs            |
|                          | Small Adult (4#) | 5401340     | 22mm OD,skid resistance,inflatable  | 50pcs            |
|                          | Medium Adult(5#) | 5401350     | 22mm OD,skid resistance, inflatable | 50pcs            |
|                          | Large Adult(6)   | 5401360     | 22mm OD,skid resistance, inflatable | 50pcs            |

Contact us for more specifications

### Type A

- ▶ Soft and flexible–transparent air cushion
- ▶ Standard connector
- ▶ Removable hook ring



| SPECIFICATION            | MODEL            | CODE NUMBER | DESCRIPTION         | Packing (carton) |
|--------------------------|------------------|-------------|---------------------|------------------|
| Anesthesia Mask (Type A) | Neonate (0#)     | 5401100     | 15mm OD,inflatable  | 30pcs            |
|                          | Infant (1#)      | 5401110     | 15mm OD, inflatable | 30pcs            |
|                          | Child (2#)       | 5401120     | 22mm OD, inflatable | 70pcs            |
|                          | Small Adult (3#) | 5401130     | 22mm OD, inflatable | 50pcs            |
|                          | Medium Adult(4#) | 5401140     | 22mm OD, inflatable | 50pcs            |
|                          | Large Adult(5#)  | 5401150     | 22mm OD, inflatable | 50pcs            |

Contact us for more specifications

Ph: 505-550-6501 & 505-565-5102, Fx: 505-814-5778

Em: sales@agsmedical.com, Web: www.agsmedical.com