



**AGS-Medical**  
**Ph: 505-550-6501 & 505-565-5102**  
**Fx: 505-814-5778**

**Em: sales@agsmedical.com, Web: www.agsmedical.com**

*Shockwave Therapy Devices List*

No.	Product Name	Model	Picture	Qualification	Quotation
1	Shock Wave Therapy Device	OCE-ESWT-001		<ul style="list-style-type: none"> <li>● CFDA</li> <li>● CE</li> <li>● FDA</li> <li>● TGA</li> <li>● ANVISA</li> </ul>	Annex 1
2	Shock Wave Therapy Device	OCE-ESWT-002		<ul style="list-style-type: none"> <li>● CFDA</li> <li>● CE</li> <li>● FDA</li> <li>● TGA</li> <li>● ANVISA</li> </ul>	Annex 2
3	Shock Wave Therapy Device	OCE-ESWT-003		<ul style="list-style-type: none"> <li>● CE</li> <li>● FDA</li> <li>● TGA</li> </ul>	Annex 3
4	Deep Muscle Stimulator	OCE-DMST-001		<ul style="list-style-type: none"> <li>● CE</li> </ul>	Annex 4
5	Ultrasonic physiotherapy equipment	OCE-ETUS-001		<ul style="list-style-type: none"> <li>● CE</li> </ul>	Annex 5





### Annex 1: Shock Wave Therapy Device (OCE-ESWT-001)

Product name	Model	Picture	Unit Price	MOQ
Shock Wave Therapy Device	OCE-ESWT-001			Q < 5
				5 ≤ Q < 10
				10 ≤ Q < 20
				Q ≥ 20
<b>Clinical Scope</b>				
<ul style="list-style-type: none"> <li>●Physiotherapy   ●Aesthetics for cellulite   ●ED treatment</li> </ul>				
<b>Technical Specifications</b>			<b>Packing</b>	
<ul style="list-style-type: none"> <li>●<b>Power Levels / Energy: 10 to 185 mJ</b> (equivalent to 1-5 bar)</li> <li>●<b>Frequency:</b> 1-22 Hz</li> <li>●<b>Modes:</b> 4 Burst modes (continuous / 4 / 8 / 12 shocks)</li> <li>●<b>Programmes:</b> 7 preset programmes, adjustable,M</li> <li>●<b>Protocols:</b> more than 25 illustrated preset treatment recommendations</li> <li>●<b>Controls:</b> Colour touch-screen for all software operations</li> <li>●<b>Applicators:</b> 6/15/25mm</li> <li>●<b>Life time:</b> Minimum <b>300,000</b> shots per applicator head; <b>4,000,000</b> shots per handpiece</li> <li>●<b>Dimensions:</b> 290×240×130mm(L/W/H)</li> <li>●<b>Weight:</b> 2.07kg (controller unit only)</li> <li>●<b>Package:</b> Aluminum Carrying Case</li> <li>●<b>Warranty:</b> 3 years for the host, 1 year for the accessories</li> </ul>				
<b>Remarks</b>				
<ul style="list-style-type: none"> <li>●MOQ for a single order quantity</li> <li>●Those who purchase shock wave therapy device from our company up to 40 sets, will get one set free from us; those who purchase accumulative total of 100 sets, will get three sets free.</li> <li>●Customers providing us information on clinical effects of our products and letters of recommendation will get a discount of US\$100 from our quotation.</li> </ul>				






## Annex 2: Shock Wave Therapy Device (OCE-ESWT-002)

Product name	Model	Picture	Unit Price	MOQ
Shock Wave Therapy Device	OCE-ESWT-002			Q < 5
				5 ≤ Q < 10
				10 ≤ Q < 20
				Q ≥ 20
<b>Clinical Scope</b>				
<ul style="list-style-type: none"> <li>●Physiotherapy</li> <li>●Aesthetics for cellulite</li> <li>●ED treatment</li> </ul>				
<b>Technical Specifications</b>			<b>Packing</b>	
<ul style="list-style-type: none"> <li>●<b>Power Levels / Energy:</b> 10 to 185 mJ (equivalent to 1-5 bar)</li> <li>●<b>Frequency:</b> 1-22 Hz</li> <li>●<b>Modes:</b> 4 Burst modes (continuous / 4 / 8 / 12 shocks)</li> <li>●<b>Programmes:</b> 7 preset programmes, adjustable,M</li> <li>●<b>Protocols:</b>more than 25 illustrated preset treatment recommendations</li> <li>●<b>Controls:</b> Colour touch-screen for all software operations</li> <li>●<b>Applicators:</b> 6/15/25mm</li> <li>●<b>Life time:</b> Minimum 300,000 shots per applicator head; 4,000,000 shots per handpiece</li> <li>●<b>Dimensions:</b>290×240×130mm(L/W/H)</li> <li>●<b>Weight:</b> 2.07kg (controller unit only)</li> <li>●<b>Package:</b> Aluminum Carrying Case</li> <li>●<b>Warranty:</b> 3 years for the host, 1 year for the accessories</li> </ul>				
<p><b>Remarks</b></p> <ul style="list-style-type: none"> <li>●MOQ for a single order quantity</li> <li>●Those who purchase shock wave therapy device from our company up to 40 sets, will get one set free from us; those who purchase accumulative total of 100 sets, will get three sets free.</li> <li>●Customers providing us information on clinical effects of our products and letters of recommendation will get a discount of US\$100 from our quotation.</li> </ul>				




### Annex 3: Shock Wave Therapy Device (OCE-ESWT-003)

Product name	Model	Picture	Unit Price	MOQ
Shock Wave Therapy Device	OCE-ESWT-003			Q < 5
				5 ≤ Q < 10
				10 ≤ Q < 20
				Q ≥ 20
<b>Clinical Scope</b>				
<ul style="list-style-type: none"> <li>●Physiotherapy   ●Aesthetics for cellulite   ●ED treatment</li> </ul>				
<b>Technical Specifications</b>			<b>Packing</b>	
<ul style="list-style-type: none"> <li>●Power Levels / Energy: <b>10 to 185 mJ</b> (equivalent to 1-5 bar)</li> <li>●Frequency: 1-22 Hz</li> <li>●Modes: continuous</li> <li>●Display mode: OLED</li> <li>●Controls: operate on handle and via bluetooth</li> <li>●Applicators: 6/15/25mm</li> <li>●Life time: Minimum <b>300,000</b> shots per applicator head; <b>4,000,000</b> shots per handpiece</li> <li>●Dimensions: 300×240×180mm(L/W/H)</li> <li>●Gross weight: 7kg</li> <li>●Package: Aluminum Carrying Case</li> <li>●Warranty: 3 years for the host, 1 year for the accessories</li> </ul>			 	
<b>Remarks</b> <ul style="list-style-type: none"> <li>●MOQ for a single order quantity</li> <li>●Those who purchase shock wave therapy device from our company up to 40 sets, will get one set free from us; those who purchase accumulative total of 100 sets, will get three sets free.</li> <li>●Customers providing us information on clinical effects of our products and letters of recommendation will get a discount of US\$100 from our quotation.</li> </ul>				




#### Annex 4: Deep Muscle Stimulator (OCE-DMST-001)

Product name	Model	Picture	Unit Price	MOQ
Deep Muscle Stimulator	OCE-DMST-001			Q < 5
				5 ≤ Q < 10
				10 ≤ Q < 20
				Q ≥ 20
<b>Clinical Scope</b>				
<b>●Physiotherapy</b>				
<b>Technical Specifications</b>			<b>Packing</b>	
<ul style="list-style-type: none"> <li>●<b>RPM:</b>1200 to 3600 rpm and adjustable</li> <li>●<b>Frequency:</b>20-60Hz</li> <li>●<b>Treatment depth:</b>20-40mm</li> <li>●<b>Applicators:</b>6/15/25mm</li> <li>●<b>Telescopic distance of head:</b>6mm</li> <li>●<b>Diameter of handle:</b>42mm</li> <li>●<b>Material of head:</b>titanium alloy</li> <li>●<b>Weight of device:</b>2.5kg</li> <li>●<b>Gross weight:</b>8kg</li> <li>●<b>Size of device:</b>240×50×136mm(L/W/H)</li> <li>●<b>Warranty:</b>1 years for controller unit</li> </ul>				
<b>Remarks</b> <ul style="list-style-type: none"> <li>●MOQ for a single order quantity</li> <li>●Those who purchase Deep Muscle Stimulator from our company up to 40 sets, will get one set free from us; those who purchase accumulative total of 100 sets, will get three sets free.</li> <li>● Customers providing us information on clinical effects of our products and letters of recommendation will get a discount of US\$ 50 from our quotation.</li> </ul>				



**Annex 5: Ultrasonic physiotherapy equipment (OCE-ETUS-001)**

Product name	Model	Picture	Unit Price	MOQ
Ultrasonic physiotherapy equipment	OCE-ETUS-001			Q < 30
				30 ≤ Q < 60
				60 ≤ Q < 100
				Q ≥ 100
<b>Clinical Scope</b>				
●Physiotherapy ●Beauty				
<b>Technical Specifications</b>			<b>Packing</b>	
<ul style="list-style-type: none"> <li>●Frequency: 1MHz / 3MHz</li> <li>●Max Intensity: 3W/cm<sup>2</sup></li> <li>●Work mode: continuous / pulsed</li> <li>●Frequency: 16Hz/48Hz/100Hz</li> <li>●Pulse duty factor: 10-90%, stepping 10%</li> <li>●RBN: &lt; 8.0</li> <li>●Ultrasound head: 5cm<sup>2</sup> (1cm<sup>2</sup> / 2cm<sup>2</sup> as option)</li> <li>●Dimension: 290 × 240 × 130mm (device)</li> <li>●Gross weight: 4kg</li> <li>●Warranty: 1 years for controller unit</li> </ul>			/	
<b>Remarks</b> <ul style="list-style-type: none"> <li>●MOQ for a single order quantity</li> <li>●Those who purchase Ultrasonic physiotherapy equipment from our company up to 60 sets, will get one set free from us; those who purchase accumulative total of 100 sets, will get two sets free.</li> <li>●Customers providing us information on clinical effects of our products and letters of recommendation will get a discount of US\$30 from our quotation.</li> </ul>				