



AGS-Medical



Skin Marker Pen

- ▶ Non-toxic
- ▶ Safe
- ▶ Resistant ink



fine tip(0.5mm)

regular tip(1.0mm)



SPECIFICATION	MODEL	CODE NUMBER	DESCRIPTION	Packing (carton)
Regular tip	Without Ruler	010110	1.0mm line width,without ruler. Individual packing.	1000pcs
	With Ruler	010111	1.0mm line width,with a ruler. Individual packing.	1000pcs
Fine tip	Without Ruler	010120	0.5mm line width,without ruler. Individual packing.	1000pcs
	With Ruler	010121	0.5mm line width, with a ruler. Individual packing.	1000pcs
Dual-tip	Without Ruler	010150	One tip is fine tip, other tip is with a ruler.	1000pcs
	With Ruler	010151	One tip is fine tip, other tip is with a ruler.	1000pcs

Contact us for more specifications

Loss of Resistance Syringe

- ▶ Consistent resistance
- ▶ Low resistance
- ▶ Safety ensured



SPECIFICATION	MODEL	CODE NUMBER	DESCRIPTION	Packing (carton)
Luer slip	5ml	1503105	Volume of syringe is 5ml	1000pcs
	7ml	1503107	Volume of syringe is 7ml	1000pcs
	10ml	1503110	Volume of syringe is 10ml	1000pcs
Luer lock	5ml	1503005	Volume of syringe is 5ml	1000pcs
	7ml	1503007	Volume of syringe is 7ml	1000pcs
	10ml	1503010	Volume of syringe is 10ml	1000pcs

Contact us for more specifications

Ph: 505-550-6501 & 505-565-5102, Fx: 505-814-5778

Em: sales@agsmedical.com, Web: www.agsmedical.com