


AGS-Medical

Ph: 505-550-6501 & 505-565-5102

Fx: 505-814-5778

Em: sales@agsmedical.com, Web: www.agsmedical.com



Name	Picture	Color	Size	Description	Packing	MOQ	FOB Price for 1x20'FCL
Nitrile gloves		blue	XS/S/M/L/XL	Powder free	100pcs/box 10boxes/ctn	100ctns	

*Our suggested size ratio - 25%S, 50%M & 25%L.

Technical Specifications:

Product Name: Nitrile Powder Free Examination Gloves

Available Colour: Blue, Black, Red, White, Purple

Application: Hospital, Beauty Salon, Working Plant, Laboratory

Raw Material: 100% Nitrile (Latex Free)

Size: XS, S, M, L, XL,

FDA Approved



Grade:AQL 1.5, AQL 2.5, AQL 4.0
Type: Powder Free & Non-Sterile

Features:

Ambidextrous, Beaded Cuff, Finger-textured, Chlorinated, Powder Free.

Standards of Conformance / Certificates:

FDA: 510K approval
ASTM : D6319, D412, D573, D5151, D6124
NFPA 1999
CE: EN455-1/2/3/4, EN420
ISO : 13458, 9001

Packing

Pack 100pcs per dispenser box, 10boxes per carton.
Measurement: 32x23x24cm
Load about 1800 cartons per 20ft container.

Performance Requirements			
Characteristic	Sampling Plan	Inspection Level	AQL
Watertight Test	Multiple Normal	G1	1.5
Visual Test	Single Normal	S4	Major 2.5
			Minor 4.0
Physical Dimension	Single Normal	S2	AQL 4.0
Physical Properties	Single Reduced	S2	AQL 4.0

Particulate Residue	NA		N = 5		≤ 2 mg/glove
Physical Dimensions					
SIZE	Extra-Small	Small	Medium	Large	Extra-Large
Palm Width ± 4mm	76	86	98	107	115
Glove Length (mm)	Minimum 230				
Thickness ±0.02mm					
Palm (center of palm)	0.1				
Finger (13 ± 3mm from tip)	0.12				
Physical Properties					
Test	Before Aging			After Aging	
Tensile Strength (MPa)	Minimum 18.0			Minimum 16.0	
Elongation (%)	Minimum 650			Minimum 500	
Force at break (N)	Minimum 9N			Minimum 6N	