



317mm*127mm*45mm L*W*H

Function introduction

- Stable and long usage life Philips special UV lamps as light source
- Small size and light weight, convenient for patients to carry
- High intensity output and obvious treatment effect
- Built-in reflector increases the efficiency of exposure
- Special comb is easy for the treatment of scalp parts
- Timer function is optional (KN-4003SL)

Specification

KN-4003A/L	1×9W UVA
KN-4003B/B	1×9W-NBLVB



UV lamp parameter

PL-S9W/10/2P	350-400nm
PL-S9W/01/2P	311-312nm

1/7/08



400mm*180mm*150mm L*W*H

Function introduction

- Stable and long usage life Philips special UV lamps as light source
- Narrow band UV light with high intensity
- Multiple configuration is optional, intelligent controlling
- Small volume, light weight, convenient to carry
- Unique structure, suitable for clinical and home use
- Optional voice prompt function(KN-4006SL)

Specification

KN-4006A/L	2×9W UVA
KN-4006B/L	2×9W-NBLVB



UV lamp parameter

PL-S9W/10/2P	350-400nm
PL-S9W/01/2P	311-312nm

1/7/08