



Innovate
Designed with 4 patents

Safety
Accord with Safety standard
9342/EEC, EN60601

Technique
Manufactured from precisely designed mould,
Stamping bed platform
Panasonic® robot welding technique
11 processes epoxy coating technique

Quality Control
34 processes for quality control


Humanized Design
Visual built-in angle
indicator with SS ball

Nice Performance
Meet various requirements for medical care,
multifunctional

Ingenuity
Exquisite conception, customization
to meet customer's requirement

World-class Branded Parts
Secure® Taiwan
Baosteel®
US Motor

With a total area of more than 30,000m², we can fulfill various requirements of production. Through the intelligent management system, each process is visualized and standardized. The manufacturing base has its strict quality management system, so every production procedure is enabled to follow their standards and every produced piece is traceable.





75°±10° Back-rest	35°±10° Knee-rest	400mm Height Adjustment 480-880mm	12° Reverse Trendelenburg	12° Trendelenburg	15° Lateral Tilt Left	15° Lateral Tilt Right	weight Weighing System 8888 kg	CPR Instant Release
----------------------	----------------------	--	---------------------------------	----------------------	-----------------------------	------------------------------	--------------------------------------	------------------------

0-90° Angle Indicator	Auto-regression
--------------------------	-----------------



Grace
Intensive Care 1

D9d

Electric ICU Bed (Lateral Tilt)
SKD-B

Technical Parameters

External Length	2280mm
External Width	1090mm
Bed Height Min.	480mm
Bed Height Max.	880mm
Mattress Platform	1925×900mm
Safe Working Load	250kg

Technical Configuration

<input checked="" type="checkbox"/> Linear Motor	3pcs
<input checked="" type="checkbox"/> Telescopic Column Motor	2pcs
<input checked="" type="checkbox"/> 6" Double Side Secure® Castors, Central Locking	4pcs
<input checked="" type="checkbox"/> ABS Groove Side Rail with Integrated LCD Control Panel	1set
<input checked="" type="checkbox"/> ABS Groove Bed Ends with Integrated LCD Control Panel	1set
<input checked="" type="checkbox"/> Ø85mm Wall Bumper Wheel	4pcs
<input checked="" type="checkbox"/> Trapeze Socket	2pcs
<input checked="" type="checkbox"/> IV Pole	1set
<input checked="" type="checkbox"/> IV Pole Prevision	4pcs
<input checked="" type="checkbox"/> Drainage Hook	2pcs
<input checked="" type="checkbox"/> Central Locking Pedal	2pcs
<input checked="" type="checkbox"/> X-ray Translucent Mattress Platform	1set

Color Options





75°±10° Back-rest	35°±10° Knee-rest	250mm Height Adjustment 570-820mm	12° Reverse Trendelenburg	12° Trendelenburg	15° Lateral Tilt Left	15° Lateral Tilt Right	weight Weighing System 8000kg	CPR Instant Release
----------------------	----------------------	--	---------------------------------	----------------------	-----------------------------	------------------------------	-------------------------------------	------------------------

0-90° Angle Indicator	Auto-regression
--------------------------	-----------------



Grace Intensive Care II

D7d

Electric ICU Bed (Lateral Tilt)
SKD-B

Technical Parameters

External Length	2230mm
External Width	1030mm
Bed Height Min.	570mm
Bed Height Max.	820mm
Mattress Platform	1960x900mm
Safe Working Load	250kg

Technical Configuration

<input type="checkbox"/> Linear Motor	2pcs
<input type="checkbox"/> Telescopic Column Motor	3pcs
<input type="checkbox"/> 6" Double Side Secure® Castors, Central Locking	4pcs
<input type="checkbox"/> ABS Grace Side Rail with Integrated Control Panel	1set
<input type="checkbox"/> ABS Grace Bed Ends with Integrated LCD Control Panel	1set
<input type="checkbox"/> Trapeze Socket	2pcs
<input type="checkbox"/> IV Pole	1set
<input type="checkbox"/> IV Pole Prevision	4pcs
<input type="checkbox"/> Drainage Hook	2pcs
<input type="checkbox"/> Central Locking Pedal	1pc

Color Options





75°±10° Back-rest	35°±10° Knee-rest	290mm H-Height Adjustment 460-750mm	12° Reverse Trendelenburg	12° Trendelenburg	CPR Instant Release	0-90° Angle Indicator



Grace Care
D8d
Electric Bed
SKD-C

Technical Parameters

External Length	2280mm
External Width	1090mm
Bed Height Min.	460mm
Bed Height Max.	750mm
Matress Platform	1925×900mm
Safe Working Load	250kg

Technical Configuration

<input checked="" type="checkbox"/> Linear Motor	4pcs
<input checked="" type="checkbox"/> 6" Double Side Secure® Castors, Central Locking	4pcs
<input checked="" type="checkbox"/> ABS Castor Side Rail with Integrated Control Panel	1set
<input checked="" type="checkbox"/> ABS Grace Bed Ends with Integrated LCD Control Panel	1set
<input checked="" type="checkbox"/> Ø285mm Wall Bumper Wheel	4pcs
<input checked="" type="checkbox"/> Trapeze Socket	2pcs
<input checked="" type="checkbox"/> IV Pole	1set
<input checked="" type="checkbox"/> IV Pole Prevision	4pcs
<input checked="" type="checkbox"/> Drainage Hook	2pcs
<input checked="" type="checkbox"/> Central Locking Pedal	2pcs

Color Options





X9x Electric Bed SKD-C

Technical Parameters

External Length	2290mm
External Width	1020mm
Bed H-Height Min.	500mm
Bed H-Height Max.	1000mm
Mattress Platform	1830-2030×820mm
Safe Working Load	250kg

Technical Configuration

Linear Motor	6pcs
6" Double Side Secure® Castors, Central Locking	4pcs
ABS Side Rail with Integrated LCD Control Panel	1 set
ABS Bed Ends One Step Removable	1 set
Ø74mm Wal Bumper Wheel	4pcs
IV Pole	1 set
IV Pole Provision	4pcs
Drainage Hook	2pcs
5" Directional Wheel Pedal	2pcs
Central Locking Pedal	2pcs
X-ray Translucent Mattress Platform	1 set
Night Light	1 set
Data Storage System	1 set
Bed Push Handle for Emergency Transfer	1 set


Color Options





C8f

Electric Bed
SKD-C



Technical Parameters

External Length	2330mm
External Width	1065mm
Bed Height Min.	540mm
Bed Height Max.	830mm
Mattress Platform	1925×900mm
Safe Working Load	250kg

Technical Configuration

<input checked="" type="checkbox"/> Linear Motor	4pcs	<input checked="" type="checkbox"/> IV Pole	1set
<input checked="" type="checkbox"/> 6" Double Side Secure® Castors, Central Locking	4pcs	<input checked="" type="checkbox"/> IV Pole Prevision	4pcs
<input checked="" type="checkbox"/> ABS Side Rail Extra Length	1set	<input checked="" type="checkbox"/> Drainage Hook	2pcs
<input checked="" type="checkbox"/> ABS Bed Ends with Integrated LCD Control Panel	1set	<input checked="" type="checkbox"/> Central Locking Pedal	1pc
<input checked="" type="checkbox"/> Ø85mm Wall Bumper Wheel	4pcs		
<input checked="" type="checkbox"/> Trapeze Socket	2pcs		



C8e

Electric Bed
SKD-C



Technical Parameters

External Length	2280mm
External Width	1065mm
Bed Height Min.	460mm
Bed Height Max.	750mm
Mattress Platform	1925×900mm
Safe Working Load	250kg

Technical Configuration

<input checked="" type="checkbox"/> Linear Motor	4pcs	<input checked="" type="checkbox"/> IV Pole	1set
<input checked="" type="checkbox"/> 6" Double Side Secure® Castors, Central Locking	4pcs	<input checked="" type="checkbox"/> IV Pole Prevision	4pcs
<input checked="" type="checkbox"/> ABS Side Rail Extra Height	1set	<input checked="" type="checkbox"/> Drainage Hook	2pcs
<input checked="" type="checkbox"/> ABS Bed Ends with Integrated LCD Control Panel	1set	<input checked="" type="checkbox"/> Central Locking Pedal	2pcs
<input checked="" type="checkbox"/> Ø85mm Wall Bumper Wheel	4pcs		
<input checked="" type="checkbox"/> Trapeze Socket	2pcs		



A6w

Electric Bed
SKD-C



Technical Parameters

External Length	2160mm
External Width	950mm
Bed Height Min.	400mm
Bed Height Max.	740mm
Mattress Platform	1925×830mm
Safe Working Load	250kg

Technical Configuration

<input checked="" type="checkbox"/> Linear Motor	1pc	<input checked="" type="checkbox"/> IV Pole Prevision	4pcs
<input checked="" type="checkbox"/> Double Drive Push Rod Motor	1pc	<input checked="" type="checkbox"/> Drainage Hook	2pcs
<input checked="" type="checkbox"/> 5" Double Side Castors	4pcs		
<input checked="" type="checkbox"/> Foldable Side Rail	1set		
<input checked="" type="checkbox"/> ABS Bed Ends with Safe Lock	1set		
<input checked="" type="checkbox"/> IV Pole	1set		



B6c

Electric Bed
SKD-C



Technical Parameters

External Length	2160mm
External Width	1060mm
Bed Height Min.	435mm
Bed Height Max.	740mm
Mattress Platform	1925×900mm
Safe Working Load	250kg

Technical Configuration

<input checked="" type="checkbox"/> Linear Motor	1pc	<input checked="" type="checkbox"/> IV Pole Prevision	4pcs
<input checked="" type="checkbox"/> Double Drive Push Rod Motor	1pc	<input checked="" type="checkbox"/> Drainage Hook	2pcs
<input checked="" type="checkbox"/> 5" Double Side Castors	4pcs		
<input checked="" type="checkbox"/> ABS Side Rail	1set		
<input checked="" type="checkbox"/> Abs Bed Ends with Safe Lock	1set		
<input checked="" type="checkbox"/> IV Pole	1set		



B4e Manual Bed SK-CD1



Technical Parameters

External Length	2160mm
External Width	1060mm
Bed Height Min.	460mm
Bed Height Max.	720mm
Mattress Platform	1925×900mm
Safe Working Load	250kg

Technical Configuration

☑ Manual Cranks	4pcs	☑ Drainage Hook	2pcs
☑ 6" Double Side Secure® Castors, Central Locking	4pcs	☑ Angle Indicator	4pcs
☑ ABS Side Rail Extra Height	1 set	☑ Central Locking Pedal	2pcs
☑ ABS Bed Ends with Safe Lock	1 set		
☑ IV Pole	1 set		
☑ IV Pole Prevision	4pcs		



F4w Manual Bed SK-CD1



Technical Parameters

External Length	2170mm
External Width	960mm
Bed Height Min.	460mm
Bed Height Max.	720mm
Mattress Platform	1925×830mm
Safe Working Load	250kg

Technical Configuration

☑ Manual Cranks	4pcs	☑ Drainage Hook	2pcs
☑ 5" Double Side Castors	4pcs	☑ Angle Indicator	4pcs
☑ Foldable Side Rail	1 set		
☑ ABS Bed Ends with Safe Lock	1 set		
☑ IV Pole	1 set		
☑ IV Pole Prevision	4pcs		



Grace Care
D3d
 Manual Bed
 SK-CD1



Technical Parameters

External Length	2230mm
External Width	1065mm
Bed Height Min.	500mm
Bed Height Max.	840mm
Mattress Platform	1925×900mm
Safe Working Load	250kg

Technical Configuration

<input checked="" type="checkbox"/> Manual Cranks	3pcs	<input checked="" type="checkbox"/> IV Pole Prevision	4pcs
<input checked="" type="checkbox"/> 5" Double Side Secure® Castors, Central Locking	4pcs	<input checked="" type="checkbox"/> Drainage Hook	2pcs
<input checked="" type="checkbox"/> ABS Grace Side Rail Extra Height	1set		
<input checked="" type="checkbox"/> ABS Grace Bed Ends with Safe Lock	1set		
<input checked="" type="checkbox"/> Central Locking Pedal	2pcs		
<input checked="" type="checkbox"/> IV Pole	1set		



Grace Care
D3w
 Manual Bed
 SK-CD1



Technical Parameters

External Length	2230mm
External Width	1015mm
Bed Height Min.	480mm
Bed Height Max.	820mm
Mattress Platform	1925×900mm
Safe Working Load	250kg

Technical Configuration

<input checked="" type="checkbox"/> Manual Cranks	3pcs	<input checked="" type="checkbox"/> Drainage Hook	2pcs
<input checked="" type="checkbox"/> 5" Covered Castors	4pcs		
<input checked="" type="checkbox"/> Foldable Side Rail	1set		
<input checked="" type="checkbox"/> ABS Grace Bed Ends with Safe Lock	1set		
<input checked="" type="checkbox"/> IV Pole	1set		
<input checked="" type="checkbox"/> IV Pole Prevision	4pcs		



A2w

Manual Bed
SK-C1



Technical Parameters

External Length	2160mm
External Width	950mm
Bed Height	500mm
Mattress Platform	1925×830mm
Safe Working Load	250kg

Technical Configuration

<input checked="" type="checkbox"/> Manual Cranks	2pcs	<input checked="" type="checkbox"/> IV Pole Provision	4pcs
<input checked="" type="checkbox"/> 5" Covered Castors	4pcs	<input checked="" type="checkbox"/> Drainage Hook	2pcs
<input checked="" type="checkbox"/> Foldable Side Rail	1set	<input checked="" type="checkbox"/> Shoes Holder	1set
<input checked="" type="checkbox"/> ABS Bed End with Safe Lock	1set		
<input checked="" type="checkbox"/> IV Pole	1set		



A3w

Manual Bed
SK-CD1



Technical Parameters

External Length	2160mm
External Width	950mm
Bed Height Min.	420mm
Bed Height Max.	740mm
Mattress Platform	1925×830mm
Safe Working Load	250kg

Technical Configuration

<input checked="" type="checkbox"/> Manual Cranks	3pcs	<input checked="" type="checkbox"/> Drainage Hook	2pcs
<input checked="" type="checkbox"/> 5" Covered Castors	4pcs		
<input checked="" type="checkbox"/> Foldable Side Rail	1set		
<input checked="" type="checkbox"/> ABS Bed End with Safe Lock	1set		
<input checked="" type="checkbox"/> IV Pole	1set		
<input checked="" type="checkbox"/> IV Pole Provision	4pcs		



T2k Manual Bed SK-C1



Technical Parameters

External Length	2130mm
External Width	970mm
Bed Height	500mm
Mattress Platform	1925×830mm
Safe Working Load	250kg

Technical Configuration

<input checked="" type="checkbox"/> Manual Cranks	2pcs	<input checked="" type="checkbox"/> IV Pole Provision	4pcs
<input checked="" type="checkbox"/> 5" Double Side Castors	4pcs	<input checked="" type="checkbox"/> Drainage Hook	2pcs
<input checked="" type="checkbox"/> Aluminium Alloy Foldable Side Rail	1set	<input checked="" type="checkbox"/> Shoes Holder	1set
<input checked="" type="checkbox"/> Aluminium Alloy Bed End	1set		
<input checked="" type="checkbox"/> IV Pole	1set		



R00 Flat Bed SK-A1

Technical Parameters

External Length	2030mm
External Width	910mm
Bed Height	500mm
Mattress Platform	1950×900mm
Safe Working Load	200kg

Technical Configuration

<input checked="" type="checkbox"/> Epoxy Coated Steel Bed Ends	2pcs
<input checked="" type="checkbox"/> IV Pole	1pc
<input checked="" type="checkbox"/> IV Pole Provision	4pcs
<input checked="" type="checkbox"/> Drainage Hook	2pcs



R10 Manual Bed Single Crank SK-B1



Technical Parameters

External Length	2030mm
External Width	910mm
Bed Height	500mm
Mattress Platform	1910×820mm
Safe Working Load	200kg

Technical Configuration

<input checked="" type="checkbox"/> Stainless Steel Bed Ends	2pcs
<input checked="" type="checkbox"/> IV Pole	1pc
<input checked="" type="checkbox"/> IV Pole Provision	4pcs
<input checked="" type="checkbox"/> Drainage Hook	2pcs
<input checked="" type="checkbox"/> Manual Cranks	1pc





H6k

Electric Homecare Bed
SKD-C



Technical Parameters

External Length	2080mm
External Width	1000mm
Bed Height Min.	460mm
Bed Height Max.	710mm
Mattress Platform	1925×830mm

Technical Configuration

<input checked="" type="checkbox"/> Linear Motor	3pcs	<input checked="" type="checkbox"/> IV Pole	1 set
<input checked="" type="checkbox"/> 6" Double Side Castors, Central Locking	4pcs	<input checked="" type="checkbox"/> IV Pole Prevision	4pcs
<input checked="" type="checkbox"/> Wood Bed Ends with Soft Leather Upholstery	1 set	<input checked="" type="checkbox"/> Drainage Hook	2pcs
<input checked="" type="checkbox"/> Wood Bed Side Panel	1 set		
<input checked="" type="checkbox"/> Aluminium Alloy Foldable Side Rail	1 set		



H2h

Manual Homecare Bed
SK-C1



Technical Parameters

External Length	2140mm
External Width	1020mm
Bed Height	400mm
Mattress Platform	1925×830mm

Technical Configuration

<input checked="" type="checkbox"/> Manual Cranks	2pcs	<input checked="" type="checkbox"/> IV Pole	1 set
<input checked="" type="checkbox"/> 5" Covered Castors	4pcs	<input checked="" type="checkbox"/> IV Pole Prevision	4pcs
<input checked="" type="checkbox"/> Wood Bed Ends	1 set	<input checked="" type="checkbox"/> Drainage Hook	2pcs
<input checked="" type="checkbox"/> Wood Bed Side Panel	1 set	<input checked="" type="checkbox"/> Shoes Holder	1pc



GB8c

Electric Orthopedic Traction Bed
SKD-C



Technical Parameters

External Length	2160mm
External Width	1060mm
Bed Height Min.	500mm
Bed Height Max.	800mm
Mattress Platform	1935×900mm

Technical Configuration

☑ Linear Motor	4pcs	☑ Drainage Hook	2pcs
☑ 5" Double Side Secure® Castors, Central Locking	4pcs	☑ Aluminium Balkan Frame	1set
☑ ABS Side Rail	1set	☑ Central Locking Pedal	1pc
☑ ABS Bed Ends with Safe Lock	1set		
☑ IV Pole Prevision	4pcs		

GS4k

Manual Orthopedic Traction Bed
SKG-4B



Technical Parameters


External Length	2130mm
External Width	960mm
Bed Height	500mm
Mattress Platform (X-ray Translucent)	1925×830mm

Technical Configuration

☑ Manual Cranks	4pcs	☑ IV Pole Prevision	4pcs
☑ 6" Double Side Castors, Central Locking	4pcs	☑ Drainage Hook	2pcs
☑ Aluminium Alloy Foldable Side Rail	1set	☑ Stainless Steel Trapeze Handle	2pcs
☑ Stainless Steel Bed Ends	1set	☑ Stainless Steel Foot Traction Frame	1set
☑ IV Pole	1pc	☑ Central Locking Pedal	1pc

GB4e

Manual Orthopedic Traction Bed
SKG-4C



Technical Parameters

External Length	2270mm
External Width	1060mm
Bed Height	500mm
Mattress Platform	1925×900mm

Technical Configuration

☑ Manual Cranks	4pcs	☑ IV Pole	1set
☑ ABS Bed Ends with Safe Lock	1set	☑ IV Pole Prevision	4pcs
☑ ABS Bed Ends with Safe Lock	1set	☑ Drainage Hook	2pcs
☑ 5" Single Side Castors, Central Locking	4pcs	☑ Stainless Steel Balkan Frame	1set
☑ Central Locking Pedal	1pc	☑ Folding Dining Table with Buffer Folding	1pc



CX2x
Children Bed
SK-C2



Technical Parameters

External Length	1800mm
External Width	900mm
Bed Height	650mm
External Height	1260mm
Mattress Platform	1760×800mm

Technical Configuration

<input checked="" type="checkbox"/> Manual Cranks	2pcs	<input checked="" type="checkbox"/> IV Pole	1pc
<input checked="" type="checkbox"/> 5" Covered Castors	4pcs	<input checked="" type="checkbox"/> Drainage Hook	2pcs
<input checked="" type="checkbox"/> Full Length Side Rail	1set	<input checked="" type="checkbox"/> Buffer Air Spring	4pcs
<input checked="" type="checkbox"/> Plastic Bed Ends with Cartoon Bear	1set		
<input checked="" type="checkbox"/> IV Pole Prevision	4pcs		

CR0q
Children Bed
SK-A2



Technical Parameters

External Length	1600mm
External Width	740mm
Bed Height	500mm
External Height	1040mm
Mattress Platform	1580×660mm

Technical Configuration

<input checked="" type="checkbox"/> 5" Covered Castors	4pcs	<input checked="" type="checkbox"/> Drainage Hook	2pcs
<input checked="" type="checkbox"/> Full Length Side Rail	1set		
<input checked="" type="checkbox"/> Integrated Bed Ends	1set		
<input checked="" type="checkbox"/> IV Pole	1pc		
<input checked="" type="checkbox"/> IV Pole Prevision	4pcs		

CT2k
Children Bed
SK-C2



Technical Parameters

External Length	1800mm
External Width	730mm
Bed Height	500mm
Mattress Platform	1580×660mm

Technical Configuration

<input checked="" type="checkbox"/> Manual Cranks	2set	<input checked="" type="checkbox"/> IV Pole	1set
<input checked="" type="checkbox"/> Aluminium Side Rail with Extra Protection	1set	<input checked="" type="checkbox"/> IV Pole Prevision	4pcs
<input checked="" type="checkbox"/> Aluminium Bed Ends with Cartoon Sticker	1set	<input checked="" type="checkbox"/> Drainage Hook	4pcs
<input checked="" type="checkbox"/> Anti-slip Bed Feet Cover	4pcs	<input checked="" type="checkbox"/> Shoes Holder	1pc



X01-1

Baby Crib
SK-AD3



Technical Parameters

External Length	875mm
External Width	525mm
Bed Height Min.	650mm
Bed Height Max.	900mm

Technical Configuration

<input checked="" type="checkbox"/> Hydraulic Pump	1 pc
<input checked="" type="checkbox"/> Bassinet	1 pc
<input checked="" type="checkbox"/> 3" Double Side Castors	4 pcs

X01-5

Baby Crib
SK-B3



X01
Stainless Steel

Powder Coated Steel

Technical Parameters

External Length	805mm±10mm
External Width	525mm±10mm
Bed Height	970mm±10mm

Technical Configuration

<input checked="" type="checkbox"/> Bassinet	1pc	<input checked="" type="checkbox"/> IV Pole	1set
<input checked="" type="checkbox"/> 3" Covered Castors	4pcs	<input checked="" type="checkbox"/> Utility Basket	1pc
<input checked="" type="checkbox"/> IV Pole Prevision	2pcs	<input checked="" type="checkbox"/> Adjustable Air Spring	1pc

Z07

Electric Bath Bed
SKD-A



Technical Parameters

External Length	1900mm
External Width	710mm
Bed Height Min.	500mm
Bed Height Max.	900mm

Technical Configuration

<input checked="" type="checkbox"/> Vertical column motor	2pcs
<input checked="" type="checkbox"/> 5" Covered Castors	4pcs
<input checked="" type="checkbox"/> Retractable basket	2pcs

Z0n

Electric Bath Bed
SKD-A



Z0y
Hydraulic

Technical Parameters

External Length	1800mm
External Width	650mm
Bed Height Min.	500mm
Bed Height Max.	730mm


Technical Configuration

<input checked="" type="checkbox"/> Linear Motor	1pc
<input checked="" type="checkbox"/> 5" Covered Castors	4pcs
<input checked="" type="checkbox"/> Stainless Steel Side Rail	1set
<input checked="" type="checkbox"/> Stainless Steel Bed Ends	1set



X12

Electric Examination Table SKD-A



Technical Parameters

External Length	1900mm
External Width	600mm
Table Height Min.	500mm
Table Height Max.	800mm

Technical Configuration

Linear Motor	1pc
Anti-slip Bed Feet Cover	4pcs
Removable Pillow	1pc

X12-1

Hydraulic Examination Table SK-BD1



Technical Parameters

External Length	1900mm
External Width	600mm
Table Height Min.	500mm
Table Height Max.	900mm

Technical Configuration

Hydraulic Pump	1pc
Adjustable Air Spring	2pcs
5" Covered Castors	4pcs
Hydraulic pedal	2sets

X26

Electric Examination Table SKD-A



Technical Parameters

External Length	1810mm
External Width	680mm
Table Height Min.	570mm
Table Height Max.	860mm

Technical Configuration

Pillow	1pc	Hand Controller	1pc
Motor	1pc	Paper Roll Holder	1pc
Mechanical Lever to Adjust Back Rest	1pc		
2" Casters	2pcs		

X27

Electric Examination Table SKD-A



Technical Parameters

External Length	1970mm
External Width	690mm
Table Height Min.	530mm
Table Height Max.	970mm

Technical Configuration

Pillow	1pc	Hand Controller	1pc
Motor	2pcs	Paper Roll Holder	1pc
Gas spring	2pcs		
4" Casters	4pcs		



X07a

Examination Table
SK-A1



X07



Technical Parameters

External Length	1900mm
External Width	600mm
Table Height	690mm

Technical Configuration

<input checked="" type="checkbox"/> Anti-slip Bed Feet Cover	4pcs

X09

Examination Table
SK-B1



X08-1



Technical Parameters


External Length	1900mm
External Width	600mm
Table Height	680mm

Technical Configuration

<input checked="" type="checkbox"/> Anti-slip Bed Feet Cover	4pcs
<input checked="" type="checkbox"/> Mechanical Lever to Adjust Back Rest	1 set

A048

Gynecological Examination Table
SKF-B



Technical Parameters

External Length	1870mm
External Width	600mm
Table Height	680mm

Technical Configuration

<input checked="" type="checkbox"/> Drawer	4pcs	<input checked="" type="checkbox"/> Lithotomy Rod	1 set
<input checked="" type="checkbox"/> Step Stool	1pc	<input checked="" type="checkbox"/> Stainless Steel Basin	1pc
<input checked="" type="checkbox"/> Retractable Foot Board	1pc	<input checked="" type="checkbox"/> Paper Roll Holder	1pc



A046

Gynecological Table SKF-B



Technical Parameters


External Length	1900mm
External Width	600mm
Table Height	800mm

Technical Configuration

<input type="checkbox"/> Crank-Wheel	2pcs	<input type="checkbox"/> Stirrup	1set	<input type="checkbox"/> Anti-slip Bed Feet Cover	4pcs
<input type="checkbox"/> Lithotomy Rod	1set	<input type="checkbox"/> Hand Bar	1set	<input type="checkbox"/> Mechanical Lever	2pcs
<input type="checkbox"/> Shoulder Holder	1set	<input type="checkbox"/> Arm Holder	1set	<input type="checkbox"/> 3" Castors with Break	4pcs

A045

Gynecological Table SKF-A



Technical Parameters

External Length	1800mm
External Width	520mm
Table Height	740mm

Technical Configuration

<input type="checkbox"/> Lithotomy Rod	1set
<input type="checkbox"/> Hand Bar	1set
<input type="checkbox"/> Mechanical Lever	2sets

A050

Gynecological Table SKF-B



Technical Parameters

External Length	1060mm
External Width	710mm
Table Height	800mm
Foot Part Dimensions	920×680×640mm

Technical Configuration

<input type="checkbox"/> Lithotomy Rod	1set
<input type="checkbox"/> 4" Castors with Break	2pcs
<input type="checkbox"/> Anti-slip Bed Feet Cover	4pcs

A045-1

Gynecological Table SKF-A



Technical Parameters

External Length	1820mm
External Width	500mm
Table Height	800mm

Technical Configuration

<input type="checkbox"/> Lithotomy Rod	1set	<input type="checkbox"/> Mechanical Lever	2sets
<input type="checkbox"/> Hand Bar	1set	<input type="checkbox"/> Stainless Steel Basin	1pc
<input type="checkbox"/> Anti-slip Bed Feet Cover	4pcs		



A99-6

Electric Gynecological Exam Couch
SKL-B



Technical Parameters

External Length	1330mm
External Width	870mm
Table Height	1450mm

Technical Configuration

<input checked="" type="checkbox"/> Simple Side Rail	1 set
<input checked="" type="checkbox"/> Leg Holder	1 set
<input checked="" type="checkbox"/> 5" Castors 4	4pcs
<input checked="" type="checkbox"/> Foot Pedal	1 pc

A99-8

Electric Gynecological Exam Couch
SKL-B



Technical Parameters


External Length	1450mm
External Width	600mm
Table Height	580mm

Technical Configuration

<input checked="" type="checkbox"/> Hand Bars	2pcs
<input checked="" type="checkbox"/> Leg Support with Cushion	2pcs
<input checked="" type="checkbox"/> Foot Padel	1 pc
<input checked="" type="checkbox"/> Stainless Steel Sewage Basin	1 pc

A99-7

Multi-Functional Electric Obstetric Bed
SKF-C



Technical Parameters

External Length	1980mm
External Width	720mm
Table Height Min.	630mm
Table Height Max.	880mm

Technical Configuration

<input checked="" type="checkbox"/> Bed Frame	1pc	<input checked="" type="checkbox"/> Hand Bars	2pcs	<input checked="" type="checkbox"/> Caster	1pc
<input checked="" type="checkbox"/> Foot Pedal	1pc	<input checked="" type="checkbox"/> Leg Support with Cushion	2pcs		
<input checked="" type="checkbox"/> Hand Controller	1pc	<input checked="" type="checkbox"/> Arm Support with Cushion	2pcs		
<input checked="" type="checkbox"/> Anesthesia Screen	1pc	<input checked="" type="checkbox"/> Stainless Steel Sewage Basin	1pc		

A105

Multi-Functional Manual Obstetric Table
SKF-C



Technical Parameters

External Length	1900mm
External Width	600mm
Table Height Min.	700mm
Table Height Max.	900mm

Technical Configuration

<input checked="" type="checkbox"/> Bed frame	1pc	<input checked="" type="checkbox"/> Leg support with cushion	2pcs
<input checked="" type="checkbox"/> Mattress	1pc	<input checked="" type="checkbox"/> Arm support with cushion	2pcs
<input checked="" type="checkbox"/> Anesthesia screen frame	1pc	<input checked="" type="checkbox"/> Stainless steel sewage basin	1pc
<input checked="" type="checkbox"/> Hand bars	2pcs		



A98-1S

Electric Obstetric Bed
SKL-A



Technical Parameters


External Length	1930mm
External Width	940mm
Table Height Min.	620mm
Table Height Max.	860mm

Technical Configuration

<input checked="" type="checkbox"/> Lithotomy rod with stirrup	1 pair	<input checked="" type="checkbox"/> Power cord	1 pc
<input checked="" type="checkbox"/> Fith basin	1 pc	<input checked="" type="checkbox"/> CD player	1 set
<input checked="" type="checkbox"/> IV pole	1 pc	<input checked="" type="checkbox"/> ABS Guardrail	1 pair
<input checked="" type="checkbox"/> Controller	1 pc	<input checked="" type="checkbox"/> Castor	1 set

A100-4A

Electric Operating Table
SKL-A



Technical Parameters

External Length	2060mm
External Width	500mm
Table Height Min.	710mm
Table Height Max.	1010mm

Technical Configuration

<input checked="" type="checkbox"/> Bed Frame	1 pc	<input checked="" type="checkbox"/> Arm Holder with Cushion	2 pcs	<input checked="" type="checkbox"/> Kidney Bridge Handle	1 pc
<input checked="" type="checkbox"/> Mattress	1 set	<input checked="" type="checkbox"/> Leg Support with Cushion	2 pcs		
<input checked="" type="checkbox"/> Anesthesia Screen Frame	1 pc	<input checked="" type="checkbox"/> Foot Support with Cushion	2 pcs		
<input checked="" type="checkbox"/> Shoulder Holder with Cushion	2 pcs	<input checked="" type="checkbox"/> Hand Controller	1 pc		



A202

Electro-Hydraulic Operating Table
SKL-A




Technical Parameters

External Length	2000mm
External Width	500mm
Table Height Min.	750mm
Table Height Max.	1050mm

Technical Configuration

<input type="checkbox"/> Bed Frame	1pc	<input type="checkbox"/> Arm Holder with Cushion	2pcs	<input type="checkbox"/> Kidney Bridge Handle	1pc
<input type="checkbox"/> Mattress	1set	<input type="checkbox"/> Leg Support with Cushion	2pcs		
<input type="checkbox"/> Anesthesia Screen Frame	1pc	<input type="checkbox"/> Foot Support with Cushion	2pcs		
<input type="checkbox"/> Shoulder Holder with Cushion	2pcs	<input type="checkbox"/> Hand Controller	1pc		



Low-position: 750mm

High-position: 1050mm

(Optional) Head Rest

Longitudinal Adjustment: 0-295mm

Longitudinal Adjustment: 0-295mm

(Optional) Lateral Arm Support

Reversed Trendelenburg: 0-22°

Trendelenburg: 0-25°

(Optional) Waste Basin For Gynecology

Kidney Bridge Elevation: 0-80mm

Anorectum Operation Position

E.N.T Operation Position

(Optional) Orthopedics Traction Frame

Lateral Tilt Right: 0-20°

Lateral Tilt Left: 0-20°



A3001

Manual Operating Table
SKT-C



Technical Parameters

External Length	2000mm
External Width	480mm
Table Height Min.	790mm
Table Height Max.	1040mm

Technical Configuration

<input type="checkbox"/> Bed Frame	1pc	<input type="checkbox"/> Arm Holder with Cushion	2pcs
<input type="checkbox"/> Mattress	1set	<input type="checkbox"/> Leg Support with Cushion	2pcs
<input type="checkbox"/> Anesthesia Screen Frame	1pc		
<input type="checkbox"/> Shoulder Holder with Cushion	2pcs		

A205

Manual Operating Table
SKT-C



Technical Parameters

External Length	2000mm
External Width	480mm
Table Height Min.	720mm
Table Height Max.	960mm

Technical Configuration

<input type="checkbox"/> Bed Frame	1pc	<input type="checkbox"/> Arm Holder with Cushion	2pcs
<input type="checkbox"/> Mattress	1set	<input type="checkbox"/> Leg Support with Cushion	2pcs
<input type="checkbox"/> Anesthesia Screen Frame	1pc		
<input type="checkbox"/> Shoulder Holder with Cushion	2pcs		

Mattress (OEM Size)



Cair Fiber



Foam



Cair Fiber+Foam



Latex



Memory Foam



TPU Composite



Air-Through Water-Proof



Sectional Foam



SKP104-1



PVC Composite



Air-Through Water-Proof



SKP002



PU Composite



Air-Through Water-Proof



SKP008



SKP006

Air Mattress



Pump

Electrical	AC220-240V/100-120V, 50HZ/60HZ
Air output	4L/min
Pressure Range	12kPa(90mmHg)
Cycle Time	5~6min

Mattress

Arrangement	2.5" bubble pad
Material	medical grade PVC
Dimension	1900×900×64mm
Weight	2.6kg

SKP012

Air Mattress



Pump

Power supply	110V/220V/50-60Hz
Air output	6-8 liter/min
Pressure range	50mmHg - 120 mmHg
Cycle time	12min
Power consumption	7W
Fuse	1A

Mattress

Standard Size	80"×32-36"
Top cover	Nylon PU /strength
Air cells	Nylon / PVC/TPU
Bottom	tarpaulin /nylon PU
Mattress weight	5.5-6.5kg
Weight Capacity	135-145kg

SKP013

Air Mattress



Pump

Power supply	110V/220V/50-60Hz
Air output	6-8 liter/min
Pressure range	50mmHg - 120 mmHg
Cycle time	12 min
Power consumption	7W
Fuse	1A

Mattress

Standard Size	80"×32-36"
Top cover	Nylon PU /strength
Air cells	Nylon / PVC/TPU
Bottom	tarpaulin /nylon PU
Mattress weight	5.5-7kg
Weight Capacity	135-145kg

SKP015

ICU Air Mattress



Pump

Power supply	110V/220V/50-60Hz
Air output	6-8 liter/min
Pressure range	50mmHg - 120 mmHg
Cycle time	12 min
Power consumption	7W
Fuse	1A

Mattress

Standard Size	80"×32-36"
Top cover	Nylon PU / strength
Air cells	Nylon / PVC / TPU
Bottom	tarpaulin / nylon PU
Mattress weight	10kg
Weight Capacity	180kg



SKS002
ABS Bedside Table

480×480×760mm(L×W×H)



Bottle Holder



With Castors



Stainless Steel Table Top

SKS008-2

ABS Bedside Table
520×430×770mm(L×W×H)



SKS005-1

SKS005-101

Steel Bedside Table
460×470×700mm(L×W×H)



SKS005

SKS006-1

Stainless Steel Bedside Table
475×390×780mm(L×W×H)



SKS020

Solid Wood Bedside Table
500×430×650mm(L×W×H)



SKS019

SKS003

ABS Bedside Table

480×480×870mm(L×W×H)



SKS008

ABS+Steel Bedside Table
530×450×700mm(L×W×H)



SKS009

MDF Bedside Table

520×430×770mm(L×W×H)



SKS013-1

ABS+Steel Bedside Table
500×570×730mm(L×W×H)



SKS009-2

Phenolic Resin Bedside Table
530×400×745mm(L×W×H)
Dining Table Max. Height: 1010mm



SKS009-3

Phenolic Resin Bedside Table
530×475×780mm(L×W×H)



SKS009-4

Phenolic Resin Bedside Table
530×475×780mm(L×W×H)



Bedside Table



SKH046-2

Overbed Table



Technical Parameters

Length	900mm
Width	450mm
Height Min.	760mm
Height Max.	1040mm

Technical Configuration

<input checked="" type="checkbox"/> Air Spring for Height Adjustment	1pc
<input checked="" type="checkbox"/> 2" Castors	4pcs

SKH046

Overbed Table



Technical Parameters

Length	760mm
Width	380mm
Height Min.	730mm
Height Max.	1130mm

Technical Configuration

<input checked="" type="checkbox"/> Air Spring for Height Adjustment	1pc
<input checked="" type="checkbox"/> Castors	4pcs

SKH242-1

U Type Overbed Table



Technical Parameters

910×450×740-1090mm

Technical Configuration

- Air Spring for Height Adjustment
- ABS Table Board



SKH242

SKH201-4

Overbed Table



Technical Parameters

800×410×730-1130mm

Technical Configuration

- Air Spring for Height Adjustment
- MDF Table Top



SKH042

Overbed Table



Technical Parameters

800×410×730-1130mm

Technical Configuration

- Air Spring for Height Adjustment
- MDF Table Top



SKH046-1

Telescopic Dining Board



Technical Parameters

870-1020×350×28mm

Technical Configuration

- Hook





SKE001

Attendant Chair

Technical Parameters

Seat Size	640×720×940mm
Bed Size	640×1900×420mm

Technical Configuration

<input checked="" type="checkbox"/> Arm	1set
<input checked="" type="checkbox"/> 2" Castors	4pcs
<input checked="" type="checkbox"/> 3" Castors	2pcs
<input checked="" type="checkbox"/> Anti-slip Feet Cover	4pcs
<input checked="" type="checkbox"/> Pillow	1pc



SKE001-3

Luxury Attendant Chair

Technical Parameters

Seat Size	710×760×950mm
Bed Size	710×1940×420mm

Technical Configuration

<input checked="" type="checkbox"/> Wood Arm Frame	1set
<input checked="" type="checkbox"/> 2" Castors	2pcs
<input checked="" type="checkbox"/> Anti-slip Feet Cover	2pcs
<input checked="" type="checkbox"/> Pillow	1pc



SKE001-4

Luxury Attendant Sofa

Technical Parameters

Seat Size	1350×900×900mm
Bed Size	1350×1950×450mm



SKE-180A

Electric Dialysis Chair

250mm Height Adjustment	0°-70° Back Section Adjustment	-75°-20° Leg Section Adjustment
-------------------------------	--------------------------------------	---------------------------------------

Technical Parameters

Upholstery Size	1900×580mm
Arm Size	450×160mm
Height	600-850mm

Technical Configuration

<input checked="" type="checkbox"/> Hand Controller	1pc
<input checked="" type="checkbox"/> Handrail Racks	1pc
<input checked="" type="checkbox"/> Powder Line	1pc
<input checked="" type="checkbox"/> Silent Caster	1set



Electric Dialysis Chair

SKE-188

Electric Dialysis Chair

200mm Height Adjustment	0-80°±5° Back Section Adjustment	-45°-25±5° Leg Section Adjustment
-------------------------------	--	---

Technical Parameters

Upholstery Size	2000×800mm
Headrest Size	390×220×100mm
Height	520-720mm

Technical Configuration

<input checked="" type="checkbox"/> Hand Controller	1pc	<input checked="" type="checkbox"/> Over Bed Table	1pc
<input checked="" type="checkbox"/> Aluminium Handrail Rail	1set	<input checked="" type="checkbox"/> IV Pole	1pc
<input checked="" type="checkbox"/> Powder Line	1pc	<input checked="" type="checkbox"/> Centre Controlled Caster	1set
<input checked="" type="checkbox"/> Head & Foot Board	1set		





SKE087

Blood Donation Chair

Technical Parameters

Seat Size	820×850×1070mm
Bed Size	820×1600×450mm



Teacup Mat



SKE091

Blood Donation Chair

Functions

Back-rest	90°±10°
Knee-rest	45°±10°
Leg Rest Down Fold	80°±10°

Technical Parameters

Seat Size	560×1030×1200mm
Bed Size	560×1650×470mm

Technical Configuration

☑ Arm	2pcs
☑ Anti-slip Feet Cover	4pcs



SKE090

Blood Donation Chair

Technical Parameters

Seat Size	650×800×450mm
-----------	---------------

Technical Configuration

☑ Arm	2pcs
☑ Anti-slip Feet Cover	4pcs



SKE004-1

Luxury Transfusion Chair

Technical Parameters

Seat Size	670×800×1050mm
-----------	----------------

Technical Configuration

☑ Anti-slip Feet Cover	4pcs
☑ Wooden Armrest	2pcs
☑ IV Pole Previsions	2pcs
☑ IV Pole	1 set



SKE005

Multi-purpose Transfusion Chair

Technical Parameters

Seat Size	750×1000×1050mm
Bed Size	750×1550×450mm

Technical Configuration

☑ Anti-slip Feet Cover	4pcs
☑ Wooden Armrest	2pcs
☑ IV Pole Previsions	2pcs
☑ IV Pole	1 set
☑ Utility Basket	1pc



Treat-waiting Chair



SKE008-1

Material: Powder Coated
1780×740×800mm(L×W×H)



SKE008

Material: Chromed Steel
1800×740×800mm(L×W×H)



SKE006-1

Material: PVC
1655×675×870mm(L×W×H)



SKE006-2

Material: PU
1765×675×905mm(L×W×H)





SKE015-2
Nurse Stool

Chrome Plated Structure with PU Foam Seat Top
Size: 400x300x660mm



SKE016
Nurse Stool

Chrome Plated Structure with PU Foam Seat Top
Size: 400x319x630mm



SKE013
Nurse Stool

Chrome Plated Structure with PU Foam Seat Top
Size: 400x400x630mm



SKE013-2
Nurse Stool

Chrome Plated Structure with PU Foam Seat Top
Size: 400x400x1116,700mm



SKE013-3
Nurse Stool

Aluminum Foot Pedal and Nylon Base
Size: 500x400x620-820mm



SKE017
Nurse Stool

Fully Stainless Steel 304, Height-Adjustable by Screws
Size: 320x460x460mm



SKE020-1
Step Stool

Size: 405x230x990x275mm



SKE051
Dining Chair

Seat and Back Made from Plastic
Size: 460x430x780mm



SKE052
Doctor Chair

Powder Coated Steel Structure
Back and Seat Section with the Elastic Fabric
Size: 600x430x800mm



SKE053
Doctor Chair

Powder Coated Steel Structure
Backrest and Seat Padded from Fabric and Sponge
Size: 500x370x800mm



SKE072
Office Chair

Mesh Cloth with Nylon Base
Size: 520x610x450-520mm



SKE073
Office Chair

Mesh Cloth with Nylon Base
Size: 480x390x420-540mm



SKE704
Office Chair

Mesh Cloth and Sponge with Steel Pedal and Base
Size: 460x590x560-760mm



SKE705
Office Chair

Mesh Cloth with Steel Base
Size: 550x480x430-550mm



SKE065
Doctor Chair

Steel Structure
Back and Seat Made of fabric and Sponge
Size: 460x560x970mm



SKE063
Doctor Chair

Steel Structure, Made of PU Leather and High-quality Sponge
Size: 580x570x930mm



SKE053-2
Training Chair

Wiring Board (370x265mm)
Size: 360x650x750mm



SKE054
Doctor Chair

Chrome Plated Structure with High Density Foam
Size: 550x620x410-510mm



SKE054-2
Office Chair

Mesh Cloth and Sponge Back, Can Be Reclined
Size: 550x620x410-510mm



SKE710
Dining Chair

Chrome Plated Structure, Made of Fabric and High-quality Sponge
Size: 540x500x800mm



SKE708
Dining Chair

Made of Wood, PU Leather and High-quality Sponge
Size: 580x580x840mm



SKE030
Wheelchair

Powder Coated Steel Structure
Seat Dimension: 395x455mm
Size: 1020x580x800mm



SKE031
Commode Chair

Height-Adjustable, with Foot Pedal
Aluminum Folding Frame, Plastic
Size: 880x130x900mm



SKE220
Rollator

Aluminum Alloy Structure, with Basket
Size: 620x690x930mm



SKE045
Meeting Chair

Size: 565mmx980mmx700mm



SKE047
Meeting Chair

Size: 565mmx980mmx700mm



SKE048
Meeting Chair

Size: 580x710x830mm



SKE049
Meeting Chair

565mmx980mmx700mm



SKZ107
Office Desk

Size: 1400x600x750mm



SKZ110
Office Desk

Size: 1200x700x750mm



SKZ111
Executive Table

Size: 1400x700x750mm



Medical Chair



SKB040(A001)

Foldable Stretcher
SKA-1C

Material: Aluminum Alloy Pipe & Canvas Fabric



Technical Parameters

Full Size	2100×500×100mm
Folded Size	580×220×210mm
Net Weight	6kg
Load Bearing	159kg



4-fold

SKB040(A004)

Foldable Stretcher
SKA-1D

Material: Aluminum Alloy & Canvas Fabric



Technical Parameters

Full Size	1840×510×210mm
Folded Size	920×510×100mm
Net Weight	7.5kg
Load Bearing	159kg



(B) 4-fold

(A) 2-fold

SKB040(A002)

Foldable Stretcher
SKA-1B

Material: Aluminum Alloy Pipe & Canvas Fabric



Technical Parameters

Full Size	2200×530×140mm
Folded Size	1100×180×130mm
Net Weight	6kg
Load Bearing	159kg



2-fold

SKB040(A006)

Scoop Stretcher
Material: Aluminum alloy



SKX-B1

SKX-B2

Technical Parameters

Full Size	2200×420×70mm
Folded Size	1660×420×70mm
Net Weight	8.5kg
Load Bearing	159kg



SKB2C01

Helicopter Stretcher
SKU-C1

Ideal for Lifting & Transporting by Helicopter



SKB2C01
SKU-C2




Technical Parameters

Full Size	2210×630×170mm
Load Bearing	272kg

Spine Board

SKB-1A



SKB2A10
1380×450×50mm

SKB2A12
1420×450×50mm

SKB2A11
1800×460×50mm

SKB2A12
1800×450×55mm

SKB2A13
1850×410×50mm

SKB1C01-1

Foldable Transfer Stair Stretcher
SKW-B2



Technical Parameters

Full Size	850×520×1370mm
Folded Size	1100×540×180mm

SKB1C02-1

Manual Stair Stretcher for Downstairs
SKW-B2



Technical Parameters

Full Size	1250×550×1630mm
Folded Size	1250×550×330mm

SKB1C02-2

Electric Stair Stretcher
SKW-B2



Technical Parameters

Full Size	1230×500×1600mm
Folded Size	1020×550×210mm



Ambulance Stretcher Trolley

Material: Aluminum Alloy Foam Mattress: PVC Cover

SKB039(A)

SKC-2B1



Technical Parameters

Full Size	1900×550×400mm	Net Weight	25.5kg
Back-rest	70°	Load Bearing	159kg

SKB039(C)

SKC-3B4



Technical Parameters

Full Size	1900×550×920mm	Net Weight	36.8kg
Lowest Size	1900×550×260mm	Load Bearing	159kg
Back-rest	70°		

SKB039(E)

SKC-3F5



Technical Parameters

Full Size	1900×540×800mm	Net Weight	45kg
Lowest Size	1900×540×200mm	Load Bearing	159kg
Back-rest	85°		

SKB039(D)

SKC-3B5



Technical Parameters

Full Size	1950×550×870mm	Net Weight	40kg
Lowest Size	1950×550×230mm	Load Bearing	159kg
Back-rest	70°		

SKB040(A)

Stretcher Trolley
SKC-1B2

Material: Aluminum alloy canvas fabric



Technical Parameters

Full Size	1950×540×880mm	Net Weight	35kg
Lowest Size	1950×540×480mm	Load Bearing	159kg

SKB038-2

Patient Trolley
SKM-E



Technical Parameters

Lowest Size	1960×700×620mm		
Highest Size	1960×700×860mm		
Back-rest	75°		

SKB037(B)

Patient Trolley
SKM-C



Technical Parameters

Full Size	1960×600×750mm		
Back-rest	75°		

SKB037(C)

Patient Trolley
SKM-G



Technical Parameters


Lowest Size	1960×700×620mm		
Highest Size	1960×700×860mm		
Back-rest	75°		



SKB041

Operation Connecting Trolley

SKK-A



Technical Parameters

Length	3650mm
Width	640mm
Min. Height	640mm
Max. Height	970mm

Technical Configuration

<input checked="" type="checkbox"/> Crank System	2pcs	<input checked="" type="checkbox"/> 6" Castors with Central Locking	8pcs	<input checked="" type="checkbox"/> Directional Wheel Release Pedal	4pcs
<input checked="" type="checkbox"/> Back Rest Air Spring	1pc	<input checked="" type="checkbox"/> IV Pole Previsions	2pcs		
<input checked="" type="checkbox"/> ABS Side Rail	1set	<input checked="" type="checkbox"/> IV Pole	1set		
<input checked="" type="checkbox"/> Directional Wheel	2pcs	<input checked="" type="checkbox"/> Central Locking Break Pedal	4pcs		



SKB041-3

Patient Transportation Trolley

SKM-F



Technical Parameters

Length	2120mm
Width	800mm
Min. Height	630mm
Max. Height	930mm

Technical Configuration

<input checked="" type="checkbox"/> USA Hydraulic Pump	2pcs	<input checked="" type="checkbox"/> Ø74mm Wall Bumper Wheel	2pcs
<input checked="" type="checkbox"/> Crank System	1set	<input checked="" type="checkbox"/> IV Pole Previsions	2pcs
<input checked="" type="checkbox"/> Back Rest Air Spring	2pcs	<input checked="" type="checkbox"/> Mattress	1pc
<input checked="" type="checkbox"/> 8" German Castors, Central Locking	4pcs	<input checked="" type="checkbox"/> IV Pole	1set
<input checked="" type="checkbox"/> Emergency Push Handle	1set	<input checked="" type="checkbox"/> Central Locking Break Pedal	2sets
<input checked="" type="checkbox"/> Directional Wheel	1pc	<input checked="" type="checkbox"/> Directional Wheel Release Pedal	2pcs
<input checked="" type="checkbox"/> Folding Side Rail	1set	<input checked="" type="checkbox"/> X-ray Translucent Bed Platform	1set



SKB041-1

Patient Transportation Trolley
SKM-F



Technical Parameters

Length	1900mm
Width	640mm
Min. Height	550mm
Max. Height	850mm

Technical Configuration

<input checked="" type="checkbox"/> Crank System	1set	<input checked="" type="checkbox"/> Directional Wheel	1pc	<input checked="" type="checkbox"/> Central Locking Break Pedal	2pcs
<input checked="" type="checkbox"/> Back Rest Air Spring	1pc	<input checked="" type="checkbox"/> IV Pole Previsions	2pcs	<input checked="" type="checkbox"/> Transfer Mattress	1pc
<input checked="" type="checkbox"/> ABS Side Rail	1set	<input checked="" type="checkbox"/> IV Pole	1set	<input checked="" type="checkbox"/> Directional Pedal	2pcs
<input checked="" type="checkbox"/> 6" Castors with Central Locking	4pcs	<input checked="" type="checkbox"/> Oxygen Cylinder Holder	1pc		

SKB041-2

Patient Transportation Trolley
SKM-F



Technical Parameters

Length	1900mm
Width	640mm
Min. Height	630mm
Max. Height	930mm

Technical Configuration

<input checked="" type="checkbox"/> Hydraulic Pump	2pcs	<input checked="" type="checkbox"/> Directional wheel	1pc	<input checked="" type="checkbox"/> Directional Pedal	2pcs
<input checked="" type="checkbox"/> Back Rest Air Spring	1pc	<input checked="" type="checkbox"/> IV Pole Previsions	2pcs	<input checked="" type="checkbox"/> Transfer Mattress	1pc
<input checked="" type="checkbox"/> ABS Side Rail	1set	<input checked="" type="checkbox"/> IV Pole	1set		
<input checked="" type="checkbox"/> 8" Castors, Central Locking	4pcs	<input checked="" type="checkbox"/> Central Locking Break Pedal	2sets		



SKB041-6

Patient Transportation Trolley
SKM-F



Technical Parameters

Length	1930mm
Width	660mm
Min. Height	510mm
Max. Height	850mm

Technical Parameters

<input checked="" type="checkbox"/> Crank System	1set	<input checked="" type="checkbox"/> Directional Wheel	1pc	<input checked="" type="checkbox"/> Oxygen Cylinder Holder	1pc
<input checked="" type="checkbox"/> Back Rest Air Spring	1pc	<input checked="" type="checkbox"/> Directional Pedal	2pcs	<input checked="" type="checkbox"/> Central Locking Break Pedal	4pcs
<input checked="" type="checkbox"/> Plastic Side Rail	1set	<input checked="" type="checkbox"/> IV Pole Previsions	2pcs	<input checked="" type="checkbox"/> Transfer Mattress	1pc
<input checked="" type="checkbox"/> 6" Castors, Central Locking	4pcs	<input checked="" type="checkbox"/> IV Pole	1set	<input checked="" type="checkbox"/> Safe Lock for Side Rail	2pcs



SKR-AB00
New All-in-One
Computer Cart



Base	ABS Plastic Shell, High-strength Metal Frame
Function	Foot Pedal Control Height Adjustment
Castors	4" Double Side Secure® Castors
Optional	Portable Plug-in High-performance Lithium Battery

SKR-B00
Simple All-in-One
Computer Cart



Base	Can be Customized Built-in Mini Computer Host
Frame	Applied For 22 Inch LCD Display
Function	Built in Power Line
Castors	4" Double Side Secure® Castors

Registered Model: SKM-A

SKR-A01
Simple All-in-one
Computer Cart



Base	ABS Plastic Shell, Built-in Power Line
Function	Handle Control Height Adjustment
Frame	Applied For 22" All in One Computer
Castors	4" Double Side Secure® Castors

SKR-IB002
Simple Laptop Cart



Function	Handle Control Height Adjustment
Frame	Applied For ≤15" Laptop
Castors	4" Double Side Secure® Castors



SKR054-ET

Emergency Trolley
Size: 750x475 x 920mm

Technical Configuration

<input checked="" type="checkbox"/> Defibrillator shelf	1pc
<input checked="" type="checkbox"/> Sliding Side Shelf	1pc
<input checked="" type="checkbox"/> Sharp Container	1pc
<input checked="" type="checkbox"/> Dust Baskets	1pc
<input checked="" type="checkbox"/> IV pole	1pc
<input checked="" type="checkbox"/> CPR Board	1pc
<input checked="" type="checkbox"/> Oxygen Bottle Holder	1pc
<input checked="" type="checkbox"/> Power Outlet	1pc
<input checked="" type="checkbox"/> Utility Container	1pc
<input checked="" type="checkbox"/> Central Lock	1set

Registered Model: SKM-B



SKR051-ET
Emergency Trolley
Size: 750×475×920mm



SKR054-AT
Anesthesia Trolley
Size: 750×475×920mm



SKR054-IV
Clinical Trolley
Size: 750×475×920mm



SKR030
Medicine Trolley
Size: 673×480×920mm



SKR001
ABS Trolley
Size: 529×366×900mm



SKR001A-1
ABS Trolley
Size: 529×366×900mm



SKR003
ABS Trolley
Size: 529×366×900mm



SKR002
ABS Trolley
Size: 529×366×900mm



SKR004
ABS Trolley
Size: 670×475×900mm



SKR007
ABS Trolley
Size: 670×475×900mm



SKR008
ABS Trolley
Size: 670×475×900mm



SKR005
ABS Trolley
Size: 670×475×900mm



SKR010-CT4
Clinical Trolley
Size: 625×475×930mm



SKR-CT723
Clinical Trolley
Size: 625×475×930mm



SKR-IV722
Infusion Trolley
Size: 625×475×930mm



SKR-MT121
Medicine Trolley
Size: 625×475×930mm



Registered Model SKM-B



SKR-ET120
Emergency Trolley
Size: 625×475×930mm



SKR-ET121
Emergency Trolley
Size: 625×475×930mm



SKR-AT120
Anesthesia Trolley
Size: 625×475×930mm



SKR-AT121
Anesthesia Trolley
Size: 625×475×930mm



SKR-ET311
Emergency Trolley
Size: 750×475×920mm



SKR-ET312
Emergency Trolley
Size: 750×475×920mm



SKR-AT311
Anesthesia Trolley
Size: 750×475×920mm



SKR-AT312
Anesthesia Trolley
Size: 750×475×920mm



SKR-MT310
Medicine Trolley
Size: 750×475×920mm



SKR-CT311
Clinical Trolley
Size: 750×475×920mm



SKR019-ITT
IV Treatment Trolley
Size: 750×475×920mm



SKR039-MT
Medicine Trolley
Size: 625×475×920mm



SKR021-5
Case History Trolley
Size: 750×410×980mm



SKR022
Case History Trolley
Size: 450×410×950mm



SKR054-MT
Medicine Trolley
Size: 750×475×920mm



SKR054-MT850
Medicine Trolley
Size: 850×520×950mm



Registered Model SKM-B

Medical Trolley



SKR026
Clinical Trolley
Size: 500×500×880mm



SKR-CT231
Clinical Trolley
Size: 780×475×1000mm



SKR-ET223
Clinical Trolley
Size: 780×475×1000mm



SKR-ET850
Emergency Trolley
Size: 850×520×990mm



SKR024-1
Case History Trolley
Size: 750×460×930mm



SKR024
Case History Trolley
Size: 435×460×930mm



SKR-CT033
Clinical Trolley
Size: 690×450×825mm



SKR-CT320
Clinical Trolley
Size: 720×470×890mm



SKR-CT151
Medicine Trolley
Size: 670×500×870mm



SKR040
Medicine Trolley
Size: 670×500×955mm



SKR-LC712
Laundry Trolley
Size: 915×460×820mm



SKR-ET053
Emergency Trolley
Size: 760×510×960mm



SKH040
Laundry Trolley
Size: 800×600×840mm



SKH040-1
Laundry Trolley
Size: 800×600×840mm



SKH027
Laundry Trolley
Size: 900×450×1050mm



SKR-WC561
Waste Collecting Trolley
Size: 800×590×890mm



Registered Model: SKP1-B



SKH006C

Instrument Trolley (Disassembly structure)
Size: 670×450×900mm

Packing volume is greatly reduced and easy to install. No welding mark, surface is much more smooth.



SKH006
Treatment Trolley
Size: 670×450×900mm



SKH005

Fan-shaped Instrument Trolley
Size: B: 1550×460×930mm
H: 1400×460×930mm
S: 1200×460×930mm



SKH015

Treatment Trolley
Size: 510×510×940mm



SKH001
Treatment Trolley
Size: 670×450×900mm



SKH002
Dressing Trolley
Size: 670×450×900mm



SKH006-1
Treatment Trolley
Size: 670×450×900mm



SKH004
Instrument Trolley
Size: 900×550×800mm



SKH029
Dressing Trolley
Size: 670×450×1000mm



SKH018
Anaesthetic Trolley
Size: 650×450×850mm



SKH032
Infusion Treatment Trolley
Size: 670×450×1700mm



SKH012
Insulated Food Cart
Heat store type



Registered Model: SKP1-B

Medical Trolley



SKH021-2
Emergency Trolley
Size: 700×400×850mm



SKH019
Emergency Trolley
Size: 750×400×850mm



SKH028
Medicine Trolley
Size: 670×450×1000mm



SKH006-101
Treatment Trolley
Size: 670×450×900mm



SKH007-1
Stainless Steel Trolley



SKH105
Wet Linen Trolley
Size: 820×600×590mm



SKH048
Oxygen Bottle Trolley



SKH088
Swab Stand
Size: 900×500×1500mm




SKH036
Stainless Steel Washing Sink
Size: 1600(W)×630(D)×1350(H)mm



Size: 1000(W)×630(D)×1350(H)mm



Size: 2000(W)×630(D)×1350(H)mm



Size: 2400(W)×630(D)×1350(H)mm

Technical Configuration

	Faucet	2pcs
	Soap Dispenser	1pc



Medical Trolley

SKH013
Stainless Steel Trolley
for Medical Record Holder (40 shelves)
Size: 630×430×900mm



SKH014
Stainless Steel Trolley
for Medical Record Holder (20 shelves)
Size: 340×400×900mm



SKH034-1
Trash Can



(1) 18L (2) 12L

SKH034-4
Trash Can



14L

SKH041 (5)
Wash Basin



SKH038-1
Mayo Trolley
Size: 600×400×870-1100mm



SKH038-2
Mayo Trolley
Size: 500×350×870-1100mm



SKH041
I.V. Pole



(1) (2)

SK-CL001 SK-CL501
Curtain Track



SKH070
Work Table
Size: 1800×700×850mm



SKH077
Work Table
Size: 1200×800/900×1600mm



SKH049
Stainless Steel Screen
Size: 2000×1870mm



SKH078
Shelf With Four Layers
Size: 1200×600×1800mm





SKN001

Stainless Steel Tank
 Φ 170mm Φ 200mm Φ 230mm
 Φ 260mm Φ 300mm



SKN002

Stainless Steel Cream Gauze Can
 Φ 80×80mm Φ 90×90mm
 Φ 100×100mm Φ 120×120mm
 Φ 141×126mm Φ 160×154mm



SKN011

Stainless Steel Sterilization Tray
 Without holes: 150×100×45mm
 200×110×45mm
 250×160×45mm
 300×200×50mm
 With holes: 150×100×45mm
 200×110×45mm
 250×160×45mm
 300×200×50mm



SKN012

Stainless Steel Tray
 300×200×48mm 350×250×48mm
 400×300×48mm 450×300×48mm
 350×500×48mm



SKN021

Stainless Steel Needle Sterilizer
 Size: 77holes 105holes



SKN022

Stainless Steel Medicine-change Bowl
 Φ 180mm Φ 200mm Φ 220mm
 Φ 240mm Φ 260mm Φ 280mm



SKN031

Stainless Steel Kidney Dish
 Shallow: L204×130mm
 M190×116mm
 S178×104mm
 Deep: L243×125mm
 M203×102mm
 S165×80mm



SKN032

Stainless Steel Sterilized Bottle/
 Medicine Cup
 Sterilized bottles:
 Φ 50×103mm Φ 50×93mm Φ 50×83mm Φ 50×53mm
 Medicine cups: Φ 44×40mm Φ 50×50mm



SKN041

Stainless Steel Soak Forceps Bucket/
 Measuring Cup/Enema Bucket
 Measuring cup: L Φ 100×150mm
 ZL Φ 130×190mm
 Enema bucket: L Φ 100×150mm
 ZL Φ 130×190mm



SKN042

Stainless Steel Forceps Bucket
 Inverted cone: Φ 80×198mm
 Φ 64×158mm
 Φ 45×113mm
 Straight: Φ 80×200mm
 Φ 64×160mm
 Φ 45×113mm



SKH090S
 Cabinet
 Size: 900×450×750mm



SKH090M
 Cabinet
 Size: 900×450×1050mm



SKH090L
 Cabinet
 Size: 900×450×1800mm



SKH091
 Cabinet
 Size: 900×450×1800mm



SKH092
 Cabinet
 Size: 900×450×1800mm



SKH093
 Cabinet
 Size: 900×450×1800mm



SKH098-2
 Locker
 Size: 900×500×1800mm



SKH098-3
 Locker
 Size: 900×500×1800mm



SKH098-4
 Locker
 Size: 900×500×1800mm



SKH098-6
 Locker
 Size: 900×500×1800mm



SKH098-9
 Locker
 Size: 900×500×1800mm



SKH059
 Single-sided
 Medicine Shelf
 Size: 800×460×2000mm



SKH060
 Rotary
 Medicine Shelf
 Size: 920×1800mm



SKH061
 Adjustable Component
 Medicine Shelf
 Size: 1030×500×2000mm



SKH062
 Two-sided Adjustable
 Medicine Tray
 Size: 800×800×2000mm

