



SKB040(A001)

Foldable Stretcher
SKA-1C

Material: Aluminum Alloy Pipe & Canvas Fabric



Technical Parameters

Full Size	2100×500×100mm
Folded Size	580×220×210mm
Net Weight	6kg
Load Bearing	159kg



4-fold

SKB040(A004)


Foldable Stretcher
SKA-1D

Material: Aluminum Alloy & Canvas Fabric



Technical Parameters

Full Size	1840×510×210mm
Folded Size	920×510×100mm
Net Weight	7.5kg
Load Bearing	159kg



(B) 4-fold

(A) 2-fold

SKB040(A002)

Foldable Stretcher
SKA-1B

Material: Aluminum Alloy Pipe & Canvas Fabric



Technical Parameters

Full Size	2200×530×140mm
Folded Size	1100×180×130mm
Net Weight	6kg
Load Bearing	159kg




2-fold

SKB040(A006)

Scoop Stretcher

Material: Aluminum alloy



Technical Parameters

Full Size	2200×420×70mm
Folded Size	1660×420×70mm
Net Weight	8.5kg
Load Bearing	159kg



SKB2C01

Helicopter Stretcher
SKU-C1

Ideal for Lifting & Transporting by Helicopter




SKB2C01
SKU-C2

Technical Parameters

Full Size	2210×630×170mm
Load Bearing	272kg



Spine Board SKB-1A



SKB2A10
1380×350×50mm

SKB2A14
1300×350×50mm

SKB2A11
1840×460×50mm

SKB2A12
1840×460×55mm

SKB2A13
1860×410×50mm

Stretcher

SKB1C01-1

Foldable Transfer Stair Stretcher
SKW-B2



Technical Parameters

Full Size	850×520×1370mm
Folded Size	1100×540×180mm

SKB1C02-1

Manual Stair Stretcher for Downstairs
SKW-B2



Technical Parameters

Full Size	1250×550×1630mm
Folded Size	1250×550×330mm

SKB1C02-2

Electric Stair Stretcher
SKW-B2



Technical Parameters

Full Size	1230×500×1600mm
Folded Size	1020×550×210mm



Ambulance Stretcher Trolley

Material: Aluminum Alloy Foam Mattress: PVC Cover

SKB039(A)

SKC-2B1



Technical Parameters

Full Size	1900×550×400mm	Net Weight	25.5kg
Back-rest	70°	Load Bearing	159kg

SKB039(C)

SKC-3B4



Technical Parameters

Full Size	1900×550×920mm	Net Weight	36.8kg
Lowest Size	1900×550×260mm	Load Bearing	159kg
Back-rest	70°		

SKB039(E)

SKC-3F5



Technical Parameters

Full Size	1900×540×800mm	Net Weight	45kg
Lowest Size	1900×540×200mm	Load Bearing	159kg
Back-rest	85°		

SKB039(D)

SKC-3B5



Technical Parameters

Full Size	1950×550×870mm	Net Weight	40kg
Lowest Size	1950×550×230mm	Load Bearing	159kg
Back-rest	70°		

SKB040(A)

Stretcher Trolley

SKC-1B2

Material: Aluminum alloy canvas fabric



Technical Parameters

Full Size	1950×540×880mm	Net Weight	35kg
Lowest Size	1950×540×480mm	Load Bearing	159kg

SKB038-2

Patient Trolley

SKM-E



Technical Parameters

Lowest Size	1960×700×620mm		
Highest Size	1960×700×860mm		
Back-rest	75°		



SKB037(B)

Patient Trolley

SKM-C



Technical Parameters

Full Size	1960×600×750mm		
Back-rest	75°		

SKB037(C)

Patient Trolley

SKM-G



Technical Parameters


Lowest Size	1960×700×620mm		
Highest Size	1960×700×860mm		
Back-rest	75°		



SKB041

Operation Connecting Trolley

SKK-A



Technical Parameters

Length	3650mm
Width	640mm
Min. Height	640mm
Max. Height	970mm

Technical Configuration

<input type="checkbox"/> Crank System	2pcs	<input type="checkbox"/> 6" Castors with Central Locking	8pcs	<input type="checkbox"/> Directional Wheel Release Pedal	4pcs
<input type="checkbox"/> Back Rest Air Spring	1pc	<input type="checkbox"/> IV Pole Previsions	2pcs		
<input type="checkbox"/> ABS Side Rail	1set	<input type="checkbox"/> IV Pole	1set		
<input type="checkbox"/> Directional Wheel	2pcs	<input type="checkbox"/> Central Locking Break Pedal	4pcs		



SKB041-3

Patient Transportation Trolley

SKM-F



Technical Parameters

Length	2120mm
Width	800mm
Min. Height	630mm
Max. Height	930mm

Technical Configuration

<input type="checkbox"/> USA Hydraulic Pump	2pcs	<input type="checkbox"/> Ø74mm Wall Bumper Wheel	2pcs
<input type="checkbox"/> Crank System	1set	<input type="checkbox"/> IV Pole Previsions	2pcs
<input type="checkbox"/> Back Rest Air Spring	2pcs	<input type="checkbox"/> Mattress	1pc
<input type="checkbox"/> 8" German Castors, Central Locking	4pcs	<input type="checkbox"/> IV Pole	1set
<input type="checkbox"/> Emergency Push Handle	1set	<input type="checkbox"/> Central Locking Break Pedal	2sets
<input type="checkbox"/> Directional Wheel	1pc	<input type="checkbox"/> Directional Wheel Release Pedal	2pcs
<input type="checkbox"/> Folding Side Rail	1set	<input type="checkbox"/> X-Ray Translucent Bed Platform	1set



SKB041-1

Patient Transportation Trolley

SKM-F



Technical Parameters

Length	1900mm
Width	640mm
Min. Height	550mm
Max. Height	850mm

Technical Configuration

<input checked="" type="checkbox"/> Crank System	1set	<input checked="" type="checkbox"/> Directional Wheel	1pc	<input checked="" type="checkbox"/> Central Locking Break Pedal	2pcs
<input checked="" type="checkbox"/> Back Rest Air Spring	1pc	<input checked="" type="checkbox"/> IV Pole Previsions	2pcs	<input checked="" type="checkbox"/> Transfer Mattress	1pc
<input checked="" type="checkbox"/> ABS Side Rail	1set	<input checked="" type="checkbox"/> IV Pole	1set	<input checked="" type="checkbox"/> Directional Pedal	2pcs
<input checked="" type="checkbox"/> 6" Castors with Central Locking	4pcs	<input checked="" type="checkbox"/> Oxygen Cylinder Holder	1pc		

SKB041-2

Patient Transportation Trolley

SKM-F



Technical Parameters

Length	1900mm
Width	640mm
Min. Height	630mm
Max. Height	930mm

Technical Configuration


<input checked="" type="checkbox"/> Hydraulic Pump	2pcs	<input checked="" type="checkbox"/> Directional wheel	1pc	<input checked="" type="checkbox"/> Directional Pedal	2pcs
<input checked="" type="checkbox"/> Back Rest Air Spring	1pc	<input checked="" type="checkbox"/> IV Pole Previsions	2pcs	<input checked="" type="checkbox"/> Transfer Mattress	1pc
<input checked="" type="checkbox"/> ABS Side Rail	1set	<input checked="" type="checkbox"/> IV Pole	1set		
<input checked="" type="checkbox"/> 6" Castors, Central Locking	4pcs	<input checked="" type="checkbox"/> Central Locking Break Pedal	2sets		



SKB041-6

Patient Transportation Trolley

SKM-F



Technical Parameters

Length	1930mm
Width	660mm
Min. Height	510mm
Max. Height	850mm

Technical Parameters

<input checked="" type="checkbox"/> Crank System	1set	<input checked="" type="checkbox"/> Directional Wheel	1pc	<input checked="" type="checkbox"/> Oxygen Cylinder Holder	1pc
<input checked="" type="checkbox"/> Back Rest Air Spring	1pc	<input checked="" type="checkbox"/> Directional Pedal	2pcs	<input checked="" type="checkbox"/> Central Locking Break Pedal	4pcs
<input checked="" type="checkbox"/> Plastic Side Rail	1set	<input checked="" type="checkbox"/> IV Pole Previsions	2pcs	<input checked="" type="checkbox"/> Transfer Mattress	1pc
<input checked="" type="checkbox"/> 6" Castors, Central Locking	4pcs	<input checked="" type="checkbox"/> IV Pole	1set	<input checked="" type="checkbox"/> Safe Lock for Side Rail	2pcs