

# ***Model: PS800-A***

---

## ***SMOKE EVACUATOR***



*Operate in  
a forest.*

### **CONTACT INFORMATION:**

AGS-Medical, Ph: +1-505-550-6501 and +1-505-565-5102

Fax: +1-505-814-5778, Email: [sales@agsmedical.com](mailto:sales@agsmedical.com)

Skype: agstech1, Web: <http://www.agsmedical.com>







# SMOKE EVACUATOR PS800-A

The PS800-A is a small, economical system designed to meet the needs of most medical specialties. It is light enough to easily move from room to room or be placed on a table.



## Features

- Provides about 99.999% filtration
- Extremely lightweight & portable
- Easy to operate
- Economical
- Simple filter changes
- Quiet

## Technical specifications

|  |   |
|--|---|
| Product                                | Smoke evacuator   |
| Model                                  | PS800-A   |
| Outlook                                | portable  |
| Filter                                 | 4 stage ULPA filtering, VLSI level                              |
| Particle sizes                         | 0.01  |
| Efficient filtering                    | 99.99995%   |
| Motor power                            | Motor power: 1100W  |
| Static evacuation mm. H <sub>2</sub> O | 1278  |
| Footswitch control                     | Available   |
| Display                                | LED display of filter condition, flow speed, footswitch control |
| Voltage                                | 220V 50/60HZ  |
| Size (cm <sup>3</sup> )                | 18x40x30  |
| Weight (Kg)                            | 20  |

## Consumable & accessory



Hose stand



Filter



Hose

