



AGS-Medical

Ph: 505-550-6501 & 505-565-5102, Fx: 505-814-5778
Em: sales@agsmedical.com, Web://www.agsmedical.com

OR & ICU SOLUTIONS

SURGICAL LIGHTS





OL-LED70/50

LED OPERATION LAMP



TECHNICAL PARAMETERS:	
Illuminance:	160000Lux/140000Lux.
Color temperature:	4500K±500K
Color rendering index(Ra. %):	≥95
Power of bulb(W):	1WxN128/1WxN64
Spot diameter(mm):	130-350
Light focus system:	Manual
Voltage:	220V 50Hz
LED service life:	60000hours
Min.installation height:	2800mm
Feature 1:	Backlight compensation function
Feature 2:	Endoscopy light source
Feature 3:	Stepless adjusting light source



MOL-LED500


Portable LED operation light

TECHNICAL PARAMETERS:	
Illuminance:	160000Lux/140000Lux.
Color temperature:	4500K±500K
Color rendering index(Ra.%):	≥95
Power of bulb(W):	1WxN128/1WxN64
Spot diameter(mm):	130-350
Light focus system:	Manual
Voltage:	220V 50Hz
LED service life:	60000hours
Min.installation height:	2800mm
Feature 1:	Backlight compensation function
Feature 2:	Endocopy light source
Feature 3:	Stepless adjusting light source



OL-LED700

Single dome LED operation light



MEL-LED11

Portable LED examination light



PARAMETERS:	MEL-LED11
Illuminance:	10,000~120000Lux
Color temperature:	4500K±500K
Color rendering index(Ra.%):	≥90
Spot diameter(d50):	65mm
Depth of illuminance:	80cm
Power of bulb(W):	1WxN20
Rated Power:	70W
LED service life:	60000hours
Voltage:	220V 50Hz
Backup battery:	optional

MEL-LED13

Portable LED examination light



PARAMETERS:	MEL-LED13
Illuminance:	80000Lux
Color temperature:	5500K±500K
High color temperature, temperature is low	
Light brightness:	adjustable
Buld power:	LED7*3W
LED service life:	50000hours
Voltage:	AC90-240V
Backup battery:	optional

OL-70/50

Operation lamp



PARAMETERS:	Halogen 70	Halogen 50
Illuminance:	40,000~160000Lux	40,000~140000Lux
Color temperature:	4200K±500K	4200K±500K
Color rendering index(Ra. %):	=97	=97
Power of bulb(W):	150WxN1	150WxN1
Spot diameter(d50-d10) (mm):	150-280	120-210
Illumination depth (mm)	=800	=800
Brightness adjustment	eight grade dimmer system	eight grade dimmer system
Number of mirrors:	3584	2150
Light focus system:	Manual	Manual
Voltage:	220V 50Hz	220V 50Hz
Input Power:	500W	
LED service life:	=1000hours	=1000hours
Min installation height:	2800mm	2800mm
Feature 1:	M&S Lights automatic switching	



MOL-50

Portable operation lamp

HOL-12/5

Hole type operation lamp



PARAMETERS:	12	5
Lamp head diameter(mm)	900	500
Illuminance:	=120,000Lux	=50,000Lux
Color temperature:	3,500K±500K	3,500K±500K
Color rendering index(Ra.%):	=88	=88
Depth of illuminance(mm):	600	600
Size of light field (mm)	160-220	160-220
Power of bulb(W):	25	25
Rated power (W)	450	
LED service life:	=300hours	=300hours
Voltage:	220V 50Hz Or 110V 60Hz	220V 50Hz Or 110V 60Hz
Min. installation height:	2,900mm	2,900mm
Feature1:	Cold light source	Cold light source



MEL-12

Portable examination light



PARAMETERS:	Hol4
Illuminance:	=40,000Lux
Color temperature:	3,700K±500K
Color rendering index(Ra.%):	=85
Spot diameter(mm):	160~220
Depth of illuminance(mm):	900
Power of bulb(W):	25
Bulb service life(hours):	1,000
Voltage:	220V 50Hz
Backup battery:	optional

MHOL-4

Portable examination light



Product colors in the catalogue may vary from actual products due to photography and printing. Please take actual produce as standard when purchasing. Please refer to actual products for standard specifications. Products specifications subject to change without notice



AGS-Medical

Ph: 505-550-6501 & 505-565-5102, Fx: 505-814-5778
Em: sales@agsmedical.com, Web://www.agsmedical.com

